

Eppendorf Exclusive - Autumn Specials

QUALITY PRODUCTS AT GREAT VALUE

- Max. rotor capacity: 18 × 1.5/2.0mL vessels
- Max. speed: 16,873 × g (14,000rpm)
- Small footprint and low access height of 230mm

Cat No.	Description	Pack qty	Promo Price £
12813222	Centrifuge 5418 with rotor FA-45-18-11 and rotor lid	Each	1,099.00
10606820	Centrifuge 5418R, refrigerated, with rotor FA-45-18-11 and rotor lid	Each	2,599.00



- High-speed mixer (3,000rpm) for fast and efficient mixing of small volumes
- Three different tube holders (0.5mL, 1.5/2.0mL, PCR 96 plates) for high vessel flexibility
- Whisper quiet operation to benefit the work environment

Cat No.	Description	Pack qty	Promo Price £
10398952	MixMate (incl. 3 tube holders: PCR 96, 0.5mL, 1.5/2.0mL)	Each	1,299.00

- High centrifugation speed of up to 20,913 × g (14,000rpm)
- Aerosol-tight Eppendorf QuickLock® caps and lids for easy, one-hand operation
- Compact footprint saves valuable bench space
- Automatic rotor recognition and imbalance detection for maximum operational safety

Cat No.	Description	Pack qty	Promo Price £
12883242	Centrifuge 5810, 230V/50-60Hz, incl. rotor A-4-62 and 15/50mL adapters	Each	3,799.00
12803252	Centrifuge 5810, 230V/50-60Hz, incl. rotor A-4-81 and 15/50mL adapters	Each	4,299.00
13447926	Centrifuge 5810, 230 V/50-60 Hz, incl. rotor S-4-104 and 15/50mL adapters for conical tubes	Each	4,599.00
12893242	Centrifuge 5810R, 230V/50-60Hz incl. rotor A-4-62 and 15/50mL adapters	Each	5,599.00
12813252	Centrifuge 5810R, 230V/50-60Hz incl. rotor A-4-81 and 15/50mL adapters	Each	5,999.00
13415356	Centrifuge 5810R, 230 V/50- 60Hz, incl. rotor S-4-104 and 15/50 mL adapters for conical tubes	Each	6,299.00



GREAT OFFERS ON RESEARCH™ PLUS PIPETTE KITS AND QUALITY SAFE-LOCK TUBES

- Eppendorf PerfectPiston®: Ultralight and highly resistant piston system
- Positions and shapes of the control buttons/displays have been optimised for stress-free work procedures
- 4-digit volume display with magnifying display for quick and accurate volume settings

Cat No.	Description	Pack qty	Promo Price £
10128494	Research Plus 3-pack option 1, 05-10, 10-100, 100-1,000uL	Each	398.00
11578663	Research Plus 3-pack option 2, 2-20, 20-200, 100-1,000uL	Each	398.00
10408474	Research Plus 3-pack option 3, 100-1000uL, 05-5mL, 1-10mL	Each	398.00



- Hinged Safe-Lock lid prevents unintentional lid opening during incubation and storage
- Polypropylene provides high resistance to chemicals, mechanical stress and temperature extremes
- Manufactured without the use of slip agents, plasticisers, and biocides

Cat No.	Description	Pack qty	Promo Price £
10292611	Safe-Lock tubes, 0.5mL, colourless	500	28.50
10318661	Safe-Lock tubes, 0.5mL, PCR clean, colourless	500	29.85
10509691	Safe-Lock tubes, 1.5mL, colourless	1,000	29.75
10398031	Safe-Lock tubes, 1.5mL, PCR clean, colourless	1,000	29.85
10038760	Safe-Lock tubes, 2.0mL, colourless	1,000	38.50
10073321	Safe-Lock tubes, 2.0mL, PCR clean, colourless	1,000	39.85

Prices valid until 31st December 2016

eppendorf

© 2016 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

Fisher Scientific
A Thermo Fisher Scientific Brand