



**Quickly identify and tally colonies with precision using this all-inclusive package designed to streamline all colony counting procedures.**

## Fisherbrand eCount Colony Counter

Combines an electronic counter with a Sharpie™ pen for marking to prevent missing or double-counting colonies. Features count verification with beeps or LED flashes, dual counting modes (up and down), a display for total counts (up to 32 separate counts), and the ability to disable the counter for writing totals and notes.

## Fisherbrand Lightbox


The Fisherbrand Lightbox serves as a LED light source for efficient bacterial colony counting. It features a spacious viewing area (125 x 100 mm) suitable for various samples, and is powered by six 1.5V AAA batteries (included).

## Fisherbrand Magnifier

The Fisherbrand Mini-Magnifier enhances object visibility by providing 1.75X magnification. It features a swiveling, molded lens for comfortable observation.

## Fisherbrand Silicone Lab Mat

This FDA food-grade silicone labmat is suitable for various laboratory settings, including cold rooms, clean rooms, and sterile suites, where high sterility is essential. Effectively safeguards countertops from temperatures up to 200°C (392°F) and demonstrates general chemical resistance to moderate oxidizing chemicals.



**Reversible**  
Blue and Gray



Fisherbrand eCount Colony Counter (15253305)

L x W x H in	L x W x H cm
7.4 x 1.2*	17 x 3*

\*Measurements shown are with Sharpie® Pen inserted

Package Includes

eCount Colony Counter	Black Sharpie® Pen (Page 48)	Hand tally Adapter
Long-life button battery included	Weighted Stand	



Fisherbrand Lightbox & Magnifier (16360666, 12816355)

Description	L x W x H in	L x W x H cm
Mini-Light Box*	6.2 x 5.9 x 0.6	15.8 x 15 x 1.4
Mini-Magnifier	6.3 x 5.3 x 8.1	16.1 x 13.5 x 20.7



Fisherbrand Silicone Lab Mat (15837944)

Colour	L x W x H in	L x W x H cm
Grey/Blue	14 x 23.5 x 0.1	35 x 60 x 0.2



Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2023 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.