



Coveralls, hooded, KLEENGUARD® A1, PPE Category I



94

NEW

The KLEENGUARD® A10 s manufactured using breathable fabric and has elasticated cuffs, ankles and waistband making for easier donning and to remove. An elasticated hood ensures a good mask fit. Blue.

Catalogue No	Size	Pack qty
SAS-960-021A	S	50
SAS-960-031U	M	50
SAS-960-041R	L	50
SAS-960-051X	XL	50
SAS-960-061L	2XL	50
SAS-960-071Y	3XL	50



Coverall, hooded, Primeguard Standard

04

Combination of soft textile feel, breathable fabric and good air permeability, coupled with good barrier properties make this an excellent choice for clean areas, laboratories etc. Limited splash protection against non-hazardous liquids. Elasticated cuffs and ankles. CE Code Type 5.

Catalogue No	Size	Pack qty
SAS-827-041L	Medium	25
SAS-827-050K	Large	25
SAS-827-060H	XL	25
SAS-827-070E	2XL	25



Coverall, hooded, Primeguard Plus

04

Ideal for use in cleanroom/laboratories where its microporous laminate provides strength with excellent comfort. Repellent to water based liquids and low hazard chemicals, very soft texture ensuring high comfort for extended periods. Good water vapour transmission rate and excellent blood barrier.

Catalogue No	Size	Pack qty
SAS-825-040C	Medium	25
SAS-825-050W	Large	25
SAS-825-060T	XL	25
SAS-825-070Q	2XL	25

