



Personal Protective Equipment

PPE Clothing Laboratory Coats

Laboratory coats, chemical splash,
Lancer style, white

91



In many laboratory applications the wrap-over style Lancer front coat provides the answer to full frontal protection. The use of chemicals in a laboratory can prove a constant hazard to standard style open neck polyester/cotton coats. The specially designed Chem-Tech material which helps repel liquid chemicals and provides the wearer with an excellent level of protection. Material: 100% polyester with a fluorochemical liquid repellent finish. CE EN 13034

Catalogue No	Size	Chest size, inch
CS849-WH-S	S	36 to 38
CS849-WH-M	M	40 to 41
CS849-WH-L	L	42 to 44
CS849-WH-XL	XL	46 to 48
CS849-WH-2XL	2XL	50 to 52

Protective coats, limited life,
KLEENGUARD® T7,
PPE Category III



94



- Suitable for laboratories, electronic industries and food services/processing
- Collared coats with pockets, back vent and press stud fastening
- Silicone free for critical areas
- EN 1513 chemical protection to parts of the body; limited use
- EN 1149-1 anti-static

Catalogue No	Size	Chest size, inch	Pack qty
SAS-817-012C	S	36 to 38	15
SAS-817-022W	M	40 to 41	15
SAS-817-032T	L	42 to 44	15
SAS-817-042Q	XL	46 to 48	15
SAS-817-052N	2XL	50 to 52	15

Laboratory coats,
Tyvek-Pro.Tech®

94



Catalogue No	Size	Chest, inch	Pack qty
TPR250-WH-L	M	40 to 41	5
TPR250-WH-M	L	42 to 44	5
TPR250-WH-XL	XL	46 to 48	5