

# Fisherbrand Carboys

## Value your lab can rely on

- Heavy wall construction provides long lasting, rugged service
- Rectangular shape optimises bench space
- Choice of materials for general use or autoclavability

- Balanced body proportions provide stability when full, prevent accidental tipping-over, optimise vertical storage space, and allow for maximum solution recovery
- Sturdy spigot (as applicable) endures the rigours of daily lab and field use
- Polypropylene (PP) carboys are autoclavable; spigot closure included with PP carboys for use when autoclaving with liquid contained inside. Follow autoclaving instructions provided with product for proper use and performance
- Metal handle securely supports the carboy when full and allows for ergonomic transport about the lab. Handle swivels out of the way when not in use. 20L size features moulded-in hand grips on the sides for added handling ergonomics and versatility
- Wide mouth closure allows for easy filling and cleaning access
- Graduations and size identification permanently moulded into the sides for efficient use
- Good chemical resistance enables use with a wide variety of lab reagents
- Materials meet FDA requirements for food use and USP VI criteria (non-toxic)
- Carboys with spigots include a hose barb insert to accommodate ¼ to 5/16-in. ID tubing. Spigot spout without the adapter fits 5/8-in. ID tubing
- Packaged and sold individually for customer convenience



Cat. No	Material	Capacity	Spigot	Autoclavable	Pack Qty
13458029	HDPE	9L	No	No	1
13478029	HDPE	20L	No	No	1
13438039	HDPE	9L	Yes	No	1
13408039	HDPE	20L	Yes	No	1
13448029	PP	9L	No	Yes	1
13468029	PP	20L	No	Yes	1
13428039	PP	20L	Yes	Yes	1
13418039	PP	9L	Yes	Yes	1
13488029	Replacement closure				6
13448039	Replacement spigot				1

For the full range of Fisherbrand products, visit:  
[eu.fishersci.com/go/fisherbrand](http://eu.fishersci.com/go/fisherbrand)



**Terms and Conditions** Please note that all sales are subject to our standard terms and conditions. Full details can be found on the inside back cover of our catalogues or visit our website [www.eu.fishersci.com](http://www.eu.fishersci.com)

Austria: (0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20  
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478  
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00  
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303  
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

**Fisher Scientific**

A Thermo Fisher Scientific Brand

UK\_EN\_2015\_11400