



Fisherbrand™
QUALITY. RELIABILITY. VALUE.

Dry Baths and Block Heaters

Isotemp™ Digital Dry Baths and Block Heaters

Isotemp™ Digital Dry Baths/Block Heaters

Versatile digital dry baths/block heaters to meet your laboratory needs

- Digital time and temperature display
- Temperature range: ambient +5°C to +130°C
- Built-in temperature sensing probe for outstanding accuracy and control
- Precise temperature control with PID circuit
- Temperature calibration offsets the temperature to desired value
- Interchangeable aluminium alloy blocks provide versatility, easy cleaning and disinfecting
- Range of configurations hold 1, 2 or 4 interchangeable modular blocks which can hold a variety of tube sizes
- Powder-coated steel body construction ensures durability
- Applications:
 - Sample preparation
 - Enzyme preservation
 - Enzyme-substrate reactions
 - DNA amplification
 - Electrophoresis gel degeneration
 - Serum coagulation
 - Heating or thawing tubes and plates





Specifications and Ordering Information

Cat. No.	15387928	15397928	15317938
Number of Blocks	1	2	4
Temperature Range	Ambient +5°C to +130°C		
Controller Type	PID digital		
Temperature Control Accuracy	At 95°C, ≤ ±0.5°C		
Temperature Uniformity	At 37°C, ≤ ±1°C		
Timing Range	0 to 99:59 minutes or continuous		
Heating-Up Time	≤ 20 minutes (30°C to 130°C @110-120V)		≤ 25 minutes (30°C to 130°C @110-120V)
Operating Temperature	10°C ~ 30°C		
Relative Humidity Tolerance	≤ 80%		
Dimensions L x W x H	288×200×100mm	318×200×100mm	450×200×100mm
Weight	2.9kg	3.3kg	4.7kg
Electrical Requirements	100-120V~50/60Hz		
Shipment Includes	Dry Bath, both UK and EU leads, block lifter		
Certification	CE cULus RoHS		
Warranty	2 years		

Accessories

Cat. No.	Description	
15327938	46 x ø6mm aluminium block	
15337938	28 x ø10mm aluminium block	
15347938	28 x 1.5mL aluminium block	
15357938	28 x 2.0mL aluminium block	
15367938	24 x ø13mm aluminium block	
15377938	15 x ø16mm aluminium block	
15387938	12 x ø18mm aluminium block	
15397938	8 x ø20mm aluminium block	
15307948	6 x ø25mm aluminium block	

Cat. No.	Description	
15317948	40 x 0.5mL aluminium block	
15327948	1 x 96-well Elisa aluminium block	
15337948	Mixed aluminium block 32 x ø6mm + 21 x ø10mm	
15347948	Mixed aluminium block 18 x 1.5mL + 10 x 2.0mL	
15357948	Mixed aluminium block 3 x ø25mm + 12 x ø12~13mm+ 6 x ø6mm	
15377948	Mixed block 30 x 0.5mL + 20 x 0.2mL	
15387948	Flat bottom, 15 x 15mL aluminium block	
15397948	Flat bottom, 4 x 50mL aluminium block	
15307958	Conical bottom, 15 x 15mL aluminium block	
15317958	Conical bottom, 4 x 50mL aluminium block	
15337958	1 x 96-well 0.2mL non-skirted PCR plate aluminium block	
15347958	1 x 96-well 0.2mL, half or full skirted PCR plate aluminium block	