



Fisherbrand™
QUALITY. RELIABILITY. VALUE.

Centrifugation

Locking lid microcentrifuge tube

Microcentrifuge tubes with locking snap cap




Convenient, pierceable locking lids snap into place with one-handed operation

- Polypropylene (PP)
- Graduated
- Flat and pierceable cap
- Locking lid to secure the sample and avoid undesired opening
- Marking area
- Autoclavable and can withstand temperatures to -86°C
- DNase/RNase and pyrogen free
- Easy one-hand use
- Applications include:
 - Sample preparation
 - Sample storage and transportation
 - Centrifugation



Locking snap cap "hook" style



	Cat. No	Capacity, mL	Colour	Bottom shape	Max. Rcf (xg)	Diameter, mm	Length, mm	Sterile	Pack qty
1	15422545	0.5	Natural	Conical	30,000	10	30	No	500
2	15432545	1.5	Natural	Conical	26,000	13	39	No	500
3	15442545	2.0	Natural	Round	25,000	13	40	No	500

Austria: (0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

Fisher Scientific

A Thermo Fisher Scientific Brand