



Disposal bin, cigarette

Smokers Cease-Fire® is the safe and practical solution to the problem of cigarette litter. Cigarette butts are safely contained in a housing specially designed to reduce the risk of fire.

Catalogue No	Colour
SAH-130-020N	Grey
SAH-130-040H	Beige

03



Disposal bins, glass, safety box

41

Fisherbrand®

Made from sturdy corrugated cardboard and designed for separating glass waste from other types of waste. Supplied flat packed to save storage space but can be assembled very quickly. Each box has a thick polyethylene bag, which acts as a liner for holding the glass waste and any liquid spills. The unique integral lid provides added safety. When the box is full the lid is removed and a rear flap is flipped up to seal the unit ready for disposal.

Catalogue No	Description	Width, mm	Depth, mm	Height, mm	Pack qty
FB55140	Medium floor model	300	300	580	6
FB55141	Large floor model	300	300	680	6
FB55142	Small bench model	200	200	260	6



Disposal bin, glass, Danibin

03

- Orange with lid and printed label
- Easily distinguished from other waste containers
- Ideal for collection and transportation of glass and crockery
- Can be disposed of at a landfill site with other domestic refuse

Catalogue No	Alt. No	Capacity, L
SAT-642-190Q	DD820	22

d.
DANIELS

