

Pouches, sterilisation



Manufactured from high quality wet strength MG bleached kraft paper for use in porous load pre-vacuum autoclaves as well as gamma sterilisers.

Double folded bottom seal for extra integrity. Opening along width. Conforms to BS 6257.

Heat seal closure

Catalogue No	Legacy No	Code	Width, mm	Depth, mm	Height, mm	Pack qty
12768975	-	B	90	50	250	1,000
12798975	-	E	190	65	250	1,000
12708985	-	F	140	50	330	1,000
12758985	-	M34	250	100	380	500
12320479	-	M36	380	125	609	250
12778985	-	V	150	125	150	1,000

Plain top closure

Catalogue No	Legacy No	Code	Width, mm	Depth, mm	Height, mm	Pack qty
12788985	-	M4	90	50	125	2,000
12798985	-	B	90	50	250	1,000
12708995	-	C	110	30	190	1,000
12718995	-	D	140	75	250	1,000
12728995	-	F	140	50	330	1,000
12738995	-	H	180	95	380	500
12748995	-	J	100	50	150	1,000
12758995	-	L	250	100	380	500
12768995	-	M	300	75	530	250
12778995	-	T	380	125	510	250
12330479	-	M35	380	125	610	250



Pouches, Steripouches and Sterireels



Manufactured from transparent polyester/polypropylene film for strength and puncture resistance, and high quality wet strength MG kraft paper.

Opens along width. Conforms to BS 6871. Suitable for both steam and ethylene oxide sterilisation.

Flat peelable pouches

Catalogue No	Legacy No	Code	Width, mm	Height, mm	Pack qty
12798995	-	P1	76	241	1,000
12709005	-	P3	102	152	1,600
12340479	-	P4	102	254	1,000
12729005	-	P5	102	298	1,000
12779005	-	P15	254	508	500
12789005	-	P16	300	510	500

Flat peelable pouches, 200m reel

Catalogue No	Legacy No	Code	Width, mm	Pack qty
12370479	-	R2	100	1
12390479	-	R4	200	1

Gusseted peelable reels, 100m length

Catalogue No	Legacy No	Code	Width, mm	Pack qty
12799045	-	RG8	100	1

