

# Bags biohazard waste



## Bag holders, biohazard, Clavies®



Recommended for use with the clear, unprinted and standard biohazard bag ranges (see above), these biohazard bag holders ensure the bag is securely supported on a stainless steel wire frame and polypropylene bottom. The bottom also serves as a containment tray in the event of a bag leakage. The tray is orange for easy identification as a biohazard containment system and can be autoclaved at 121°C (250°F) to maintain sterility. For ease of transportation they are supplied unassembled, however assembly only takes a few minutes and once done becomes permanent.

Catalogue No	Legacy No	Tray dimensions, mm	Height, mm	Bag size, mm
11322154	-	185 x 356	470	300 x 610
11332154	-	356 x 356	546	610 x 760
11342154	-	356 x 356	696	610 x 910



## Bags, biohazard, autoclavable, HDPE and bag holder



For convenient biohazard disposal, recommended for use with the 10633305 desk size bag holder. The high molecular weight, high density polyethylene bags are bright red and printed with the biohazard symbol, with a sterilisation indicator patch that turns black at 121°C. Opening along width.

### Bags

Catalogue No	Legacy No	Width, mm	Height, mm	Thickness, µm	Pack qty
10013284	-	220	280	38	100

### Bag holder

Catalogue No	Legacy No	Description
10633305	-	Desk size rectangular epoxy coated wire holder



## Bag clips, Clavies®



Designed to securely close autoclave bags that are filled with waste material prior to being autoclaved.

Catalogue No	Legacy No	Alt. No	Description	Pack qty
11382154	-	F131900100	Polypropylene clips, orange, for 300mm x 610mm to 480mm x 580mm bags	100