

Basins, porcelain, evaporating, round bottom, with spout



Shallow form. Sizes up to 175mL glazed except rim, 1,000mL and larger are glazed except outside base.

Catalogue No	Legacy No	Capacity, mL	Diameter, mm	Height, mm
10128160	-	50	75	26
10443001	-	75	80	32
10702242	-	95	85	36
10178250	-	115	100	40
10588922	-	140	105	42
10678181	-	175	110	45
10483381	-	285	125	50
10483571	-	420	150	45
10732622	-	580	160	64
10628191	-	1,000	200	80
10628351	-	2,000	250	100

Basins, PTFE, evaporating, Fisherbrand®



Manufactured from pure PTFE with a smooth internal finish.

Inert and usable to 280°C.

Tall form with spout

Catalogue No	Legacy No	Capacity, mL	Diameter, mm	Height, mm
10196710	-	25	38	34
10577092	-	50	40	40
10577282	-	100	64	64
10034330	-	150	78	78
10287900	-	180	80	80
10646921	-	250	97	97
10517482	-	350	100	100

Flat form without spout

Catalogue No	Legacy No	Capacity, mL	Diameter, mm	Height, mm
10572202	-	25	42	25
10572392	-	50	60	20
10088171	-	100	80	28
10098171	-	180	80	46
10552972	-	350	100	55

Flat form with spout

Catalogue No	Legacy No	Capacity, mL	Diameter, mm	Height, mm
10383592	-	25	65	12
10582392	-	50	82	15
10222142	-	100	102	20
10110432	-	180	130	22
10028191	-	350	132	36
10028201	-	400	136	40

Are you online?

- See the price you pay
- Live stock availability
- Order tracking
- Fully secure online ordering



International: www.fisher.co.uk