

Beakers polypropylene



Beakers, polypropylene, squat form, ISO 7056, Thermo Scientific Nalgene®



Tapered with ribs for easy stacking. Clear graduations. Autoclavable.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Pack qty
10116310	-	1201-0030	30	12
10556361	-	1201-0050	50	12
10176350	-	1201-0100	100	12
10340481	-	1201-0150	150	12
10547882	-	1201-0250	250	6
10166400	-	1201-0400	400	6
10761411	-	1201-0600	600	4
10084810	-	1201-1000	1,000	3
10084820	-	1201-2000	2,000	1
10094820	-	1201-4000	4,000	1



Beakers, polypropylene, tapered, moulded graduations, ISO 7056



Azlon

- Tapered stacking beaker designed and manufactured to ISO 7056
- Moulded graduations and precise non-drip spout, allowing easy pouring
- Autoclavable up to 121°C

Catalogue No	Legacy No	Alt. No	Capacity, mL	Subdivisions, mL	Height, mm	O.D., mm	Pack qty
11303474	-	BPM0025P	25	5	51	41.5	10
11313474	-	BPM0050P	50	10	61	50.5	10
12505776	-	BPM0100P	100	20	73	61.5	10
11333474	-	BPM0250P	250	50	96	80.5	10
11343474	-	BPM0500P	500	50	121	102.5	10
11353474	-	BPM1000P	1,000	100	148	120	5
11363474	-	BPM2000P	2,000	200	188	163	1
11373474	-	BPM3000P	3,000	200	214	172	1
12019837	-	BPM5000P	5,000	500	250	230	1



Beakers, polypropylene, straight sided, translucent, ISO 7056



Azlon

Excellent chemical resistance, and autoclavable to 121°C.

With spout and blue, printed graduations. Not suitable for use with hotplates.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Subdivisions, mL	O.D., mm	Pack qty
10237974	-	BDA 208P	10	1	32	10
10064274	-	BDA 212P	25	2.5	42	10
10799375	-	BDA 218P	50	5	51	10
10452205	-	BDA 224P	100	10	62	10
10156444	-	BDA 230P	250	25	82	10
10330045	-	BDA 234P	400	50	94	10
10074184	-	BDA 236P	500	50	100	10
10176304	-	BDA 238P	600	50	106	5
10349324	-	BDA 240P	1,000	100	127	5
10329564	-	BDA 242P	2,000	200	158	1