

Beakers

plastic, tri-cornered



Beakers, polypropylene, tri-cornered, Fisherbrand®



Fisherbrand®

Suitable for use with commonly used acids, alkalis and solvents.

Each beaker has three drip-free pouring spouts.

Moulded graduations, stackable.

Catalogue No	Legacy No	Capacity, mL	Subdivisions, mL	Height, mm	O.D., mm	Pack qty
11749398	-	50	5	57	49	100
11759398	-	100	10	72	58	100
11769398	-	250	10	90	76	100
11779398	-	400	20	108	85	100
11789398	-	800	50	133	107	100
11799398	-	1,000	50	145	115	100



Beakers, polypropylene, Tri-Stir



Each beaker has a 'V' notch at the corners for easy pouring, easy to read graduations in ounces and millilitres.

Flat bottom surface that will keep a stirrer bar centered. Stackable. Autoclavable.

Catalogue No	Legacy No	Capacity, mL	Height, mm	I.D., mm	Pack qty
11378013	-	50	58	43 (neck) 36 (base)	10
11388013	-	100	74	53 (neck) 44 (base)	10
11398013	-	250	93	72 (neck) 60 (base)	10
11308023	-	400	112	81 (neck) 67 (base)	10
11318023	-	800	140	103 (neck) 86 (base)	10
11328023	-	1,000	150	111 (neck) 91 (base)	10



Beakers, PTFE, squat form, Fisherbrand®



Fisherbrand®

Non-contaminating and inert to highly corrosive reagents.

- Non-wettable, lightweight and heat resistant
- Can be heated on a hotplate up to 280°C
- Easy to clean
- Spout for easy pouring
- Isostatically moulded from pure PTFE

Catalogue No	Legacy No	Capacity, mL
10733271	-	5
10166570	-	10
10288380	-	25
10628461	-	50
10792921	-	100
10713101	-	250
10349740	-	500
10044630	-	1,000
10762751	-	2,000
10237860	-	5,000
10248480	-	150