

Caps, polyvinyl lined plastic, white, Fisherbrand®



General purpose suitable for a wide range of applications. Good chemical resistance suitable for mild acids, alkalis, solvents, alcohols, oils and aqueous products; poor for active hydrocarbons and bleaches.

1mm polyvinyl film bonded to 1mm HDPE on a #30 white pulp paper backing.



Superior to plain pulp paper because it provides excellent moisture barrier.

For Fisherbrand® Boston round and wide mouth packer bottles.

Catalogue No	Legacy No	Thread style	Pack qty
11583792	-	18/400	576
11593792	-	20/400	576
11503802	-	22/400	576
11513802	-	24/400	576
11523802	-	28/400	576
11533802	-	28/400	3,000
11543802	-	33/400	576
11553802	-	33/400	2,100
11563802	-	38/400	576
11583802	-	43/400	576
11593802	-	43/400	1,760
11503812	-	45/400	576
11513812	-	45/400	1,520
11523812	-	48/400	576
11533812	-	53/400	288
11543812	-	53/400	1,080
11558162	-	58/400	288
11568162	-	63/400	288
11788699	-	70/400	288
11588162	-	70/400	640
11598162	-	83/400	144
11508172	-	89/400	144

For Fisherbrand® wide mouth packer bottles, please see pages 176 to 178.



Caps, Poly-Seal® (polycone) lined plastic, black, Fisherbrand®



Unique problem-solving type of liner. This liner is stress crack resistant and offers superior torque retention and excellent sealing characteristics. Manufactured from LDPE. Unique cone design provides a wedge type seal that not only seals across the top but also across the inside diameter. However, it is recommended the liners be tested for leak seal prior to use.



For Fisherbrand® Boston round and wide mouth packer bottles.

Catalogue No	Legacy No	Thread style	Pack qty
11518172	-	18/400	576
11528172	-	20/400	576
11538172	-	20/400	5,500
11548172	-	22/400	576
11798699	-	24/400	576
11578172	-	24/400	4,200
11598172	-	28/400	3,100
11718709	-	33/400	576
11518182	-	33/400	2,300
11528182	-	38/400	576
11538182	-	38/400	1,600

Caps, PTFE faced PE lined plastic, white, Fisherbrand®



Typical applications include analytical laboratory samples, high purity chemicals, strong acids and solvents. Excellent for environmental samples, pharmaceuticals and diagnostic reagents.



PTFE faced foamed polyethylene liner that offers the excellent chemical resistance of PTFE with the compressibility and sealing properties of polyethylene foam.

For Fisherbrand® Boston round and wide packer bottles.

Catalogue No	Legacy No	Thread style	Pack qty
11548182	-	18/400	500
11558182	-	20/400	500
11568182	-	22/400	500
11578182	-	24/400	500
11588182	-	28/400	500
11598182	-	28/400	3,000
11508192	-	33/400	500
11518192	-	33/400	2,100
11528192	-	38/400	500
11538192	-	43/400	500
11548192	-	45/400	500
11558192	-	48/400	500
11568192	-	53/400	500
11578192	-	58/400	500
11588192	-	63/400	288
11598192	-	70/400	144
11508202	-	83/400	144
11518202	-	89/400	144

FINDING THE PRODUCTS YOU NEED