

Bottles

glass reagent



Bottles, borosilicate glass, laboratory, reagent and media, crimp top



Ideal for reaction vessels, sample storage and freeze drying applications.

Conforms to EP Type 1 standards.

Clear or amber bottles with crimp top and a wide range of seals and stoppers. Autoclavable.

Clear

Catalogue No	Legacy No	Alt. No	Capacity, mL	I.D. x O.D., mm	Dimensions [dia. x h], mm	Pack qty
10520853	-	223738	5	13 x 20 (neck)	23 x 47	288
10768134	-	223739	10	13 x 20 (neck)	25 x 54	288
10453712	-	223742	20	13 x 20 (neck)	32 x 58	288
10632973	-	223743	30	13 x 20 (neck)	37 x 63	288
10642973	-	223745	50	13 x 20 (neck)	43 x 73	288
10021362	-	223747	100	13 x 20 (neck)	52 x 95	144
10071382	-	223748	125	13 x 20 (neck)	54 x 107	144

Amber

Catalogue No	Legacy No	Alt. No	Capacity, mL	I.D. x O.D., mm	Dimensions [dia. x h], mm	Pack qty
10424292	-	223760	5	13 x 20 (neck)	23 x 47	288
10434292	-	223761	10	13 x 20 (neck)	25 x 54	288
10600173	-	223762	20	13 x 20 (neck)	32 x 58	288
10071392	-	223763	30	13 x 20 (neck)	27 x 63	288
10464672	-	223764	50	13 x 20 (neck)	43 x 73	288
10011422	-	223766	100	13 x 20 (neck)	52 x 95	144

Rubber stoppers - straight plug style

Catalogue No	Legacy No	Alt. No	Description	For neck dimensions, mm	Pack qty
10021492	-	224100-172	Natural red rubber	13 x 20 (I.D. x O.D.)	1,000
12811144	-	224100-173	Grey chlorobutyl/45	13 x 20 (I.D. x O.D.)	1,000
10403085	-	224100-175	PTFE coated grey chlorobutyl-isoprene blend/40	13 x 20 (I.D. x O.D.)	1,000

Rubber stoppers - 2 leg lyophilisation

Catalogue No	Legacy No	Alt. No	Description	For neck dimensions, mm	Pack qty
10452242	-	224100-192	Grey chlorobutyl-isoprene blend/50	13 x 20 (I.D. x O.D.)	1,000

Rubber stoppers - 3 leg lyophilisation

Catalogue No	Legacy No	Alt. No	Description	For neck dimensions, mm	Pack qty
11764199	-	224100-202	Grey chlorobutyl/50	13 x 20 (I.D. x O.D.)	1,000

Aluminium crimp seals

Catalogue No	Legacy No	Alt. No	Description	For neck dimensions, mm	Pack qty
10270322	-	24178-01	Open top, natural	20 (O.D.)	1,000
10194172	-	224183-01	Centre disc tear out, natural	20 (O.D.)	1,000
10721504	-	224191	Solid top, natural	20 (O.D.)	1,000
10633543	-	224193-01	Tear off, natural	20 (O.D.)	1,000