

Bottles, borosilicate glass, plastic coated, reagent and media, screw neck, ISO 4796, with retrace code



The plastic coating protects against mechanical damage and contains the contents should the glass break. The co-polymer based plastic material remains transparent and adheres to the glass surface even after repeated sterilisation in an autoclave. Steam up to 135°C; do not use distilled water for cooling. The plastic coating is largely UV absorbent to 380nm so the bottles can be used in a similar way to amber bottles. Graduated without cap or pouring ring. Do not autoclave bottles with a tightly screwed cap.

For caps and pouring rings, see page 175.

Catalogue No	Legacy No	Alt. No	Capacity, mL	ISO thread GL	Height, mm	O.D., mm	Pack qty
13158083	-	1092676	25	25	70	36	10
13151544	-	1092677	50	32	87	46	10
13127384	-	2180524	100	45	100	56	10
13157216	-	218052901	150	45	115	62	10
13128033	-	2180536	250	45	138	70	10
13137374	-	2180544	500	45	176	86	10
13167216	-	218055101	750	45	208	95	10
12387178	-	2180554	1,000	45	225	101	10
12811044	-	2180563	2,000	45	262	136	1
13177216	-	218056903	3,500	45	300	160	1
13128083	-	2180573	5,000	45	335	186	1
13187216	-	218056804	10,000	45	415	227	1
13197216	-	218058801	15,000	45	450	268	1
13107226	-	218059103	20,000	45	510	288	1



Bottles, square, borosilicate, reagent and media, screw neck, DURAN®, with retrace code



With all the properties of the DURAN® round bottles, plus the additional benefit of saving space. Using a square bottle can save up to 20.3% in space. Supplied with polypropylene caps and pouring rings.

For spare caps and pouring rings, see page 175.

Catalogue No	Legacy No	Alt. No	Capacity, mL	ISO thread GL	Height, mm	O.D., mm	Pack qty
13184594	-	21820245	100	32	109	50	10
13157614	-	21820365	250	45	143	64	10
13144584	-	21820445	500	45	181	78	10
13177134	-	21820545	1,000	45	222	94	10



Bottles, amber and clear borosilicate glass, GLS 80 wide screw neck, ISO 3685, with graduations and retrace code



Manufactured from borosilicate 3.3 glass to ISO 3685.

These bottles have a wide screw thread neck with an inside diameter of 65mm, allowing easy and safe filling, mixing and removal of powders, media, granules and fluids. Each bottle has a polypropylene GLS 80 quick release closure and is printed with a retrace code allowing web access batch certificates of conformity. All bottles are available with Duran® Protect plastic coating for additional safety. See our website for further information.

Clear

Catalogue No	Legacy No	Alt. No	Capacity, mL	ISO thread GL	Height, mm	O.D., mm	Pack qty
13127236	-	218603656	250	80	111	95	10
13164554	-	1112627	500	80	148	101	10
13127254	-	1112713	1,000	80	218	101	10
13197224	-	1112715	2,000	80	248	136	10
13137236	-	218606953	3,500	80	275	160	1
13144444	-	1113949	5,000	80	311	182	1
13131594	-	1113950	10,000	80	395	227	1
13127544	-	1113951	20,000	80	483	288	1
13147236	-	1200265	30,000	80	553	340	1
13157236	-	1200154	50,000	80	587	400	1



entry continued