

## Accessories and spare parts for GL25, GL32, GL45 thread bottles



Please see our website for further additions to the range.

Polypropylene screw caps, with lip seal, autoclavable to 140°C

Catalogue No	Legacy No	Alt. No	Colour	ISO thread GL	Pack qty
13107644	-	2923913	Blue	25	10
13177444	-	2923919	Blue	32	10
12307188	-	2923928	Blue	45	10
13137154	-	293382802	Yellow	45	10
13107236	-	293382868	Green	45	10
13117236	-	293382884	Grey	45	10



Polypropylene pouring rings, autoclavable to 140°C (for screw caps above)

Catalogue No	Legacy No	Alt. No	Colour	ISO thread GL	Pack qty
13197644	-	2924219	Blue	32	10
13138033	-	2924228	Blue	45	10
13147226	-	1089911	Green	45	10
13157226	-	1089914	Grey	45	10
13167226	-	1089917	Yellow	45	10



PBT screw caps, with PTFE protected seals, autoclavable to 180°C, red

Catalogue No	Legacy No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
13144354	-	2924013	-	25	10
13181764	-	2924019	-	32	10
13197154	-	2924028	-	45	10



EFTE pouring rings, autoclavable to 180°C, red (for screw caps above)

Catalogue No	Legacy No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
13184344	-	2924428	-	45	10

PBT caps with aperture, autoclavable to 180°C, red

Catalogue No	Legacy No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
13197434	-	2922708	20	32	10
13104264	-	2922710	34	45	10



Septa seals, silicone (for piercing) for aperture caps

Catalogue No	Legacy No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
13187164	-	2924608	29	32	100
13177374	-	2924610	42	45	100

## Premium cap, GL45 thread



Manufactured from thermoplast/duraplast (PTch260), this cap has the highest chemical resistance.

The cap contains no colourants, and has a welded-in silicone seal coated with PTFE on both sides. Matching pouring rings are available and guarantee non-dripping and safe handling. The premium cap has a temperature resistance of -196°C to 260°C and therefore suitable for many sterilisation methods such as autoclaving and dry sterilisation. The premium cap meets all the requirements of both US pharmacopoeia (USP class VI) and European pharmacopoeia (EP) and therefore complies with FDA standards. Autoclavable.



Caps

Catalogue No	Legacy No	Alt. No	ISO thread GL	Pack qty
13177226	-	1129600	25	5
13144324	-	1088679	45	5

Pouring rings

Catalogue No	Legacy No	Alt. No	ISO thread GL	Pack qty
13181864	-	1088678	45	5

For membrane screw caps, ideal for the storage or transport of gas generating media and autoclaving, please visit our website.

