

Bottles, HDPE, carboys, rectangular, Thermo Scientific Nalgene®



Space-saving design with moulded graduations.

Built-in shoulder loops have convenient stainless steel bail handle attached. 20L size has moulded-in handgrips. Polypropylene screw closure. Leakproof.

Catalogue No	Legacy No	Alt. No	Capacity, L	Height, mm	Neck I.D., mm	Closure size, mm
10706801	-	2211-0020	9	358	86	100
10180341	-	2211-0050	20	384	86	100



Bottles, fluorinated HDPE, translucent, Thermo Scientific Nalgene®



Safe, durable, low cost range with enhanced long term container performance.

Fluorinated inside and outside to reduce solvent absorption and penetration. Useful for most acids, alkalis and aggressive organic solvents. With fluorinated polypropylene screw closure. Leakproof.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10702532	-	2097-0008	250	17	133	61	24	12
10382991	-	2097-0016	500	21	171	73	24	12
10643121	-	2097-0032	1,000	27	216	92	24	6
10323141	-	2097-0005	2,000	27	246	122	24	1
10159360	-	2097-0010	4,000	27	334	153	24	1
10750982	-	2097-0020	10,000	64	389	249	24	1
10383231	-	2097-0050	20,000	64	513	282	24	1



Bottles, polypropylene, carboys, wide neck, Thermo Scientific Nalgene®



Excellent for large batches of culture media, distilled water and other solutions.

Graduated TPE gasket. Polypropylene screw closure. Leakproof.

Catalogue No	Legacy No	Alt. No	Capacity, L	Neck size, mm	Height, mm	O.D., mm	Closure size, mm
10413491	-	2250-0020	10	64	389	249	83B
10507701	-	2250-0050	20	64	513	282	83B
10343231	-	2250-0130	50	64	673	376	83B



Closures, filling/venting, Thermo Scientific Nalgene®



For aseptic transfer of media, biological reagents, pure water and chemicals from and to carboys and bottles.

Each carboy/bottle compatible with this range indicates the closure size, either 53mm or 83mm. Includes two pieces of Thermo Scientific Nalgene 860 tissue culture tubing for aspiration, dispensing and filling splash guard. Bottle not included.

Catalogue No	Legacy No	Alt. No	Closure	Height, mm	O.D., mm	I.D., mm
10589032	-	2162-0531	53B	69	67	6
10350911	-	2162-0830	83B	98	102	13
10686941	-	2162-0831	83B	98	102	6

