

Bottles, polypropylene, heavy duty, wide neck



Azlon

Clear and amber high quality, heavy duty bottles for easy filling and emptying of solid or liquid samples. Non-stick natural polypropylene screw caps. Excellent general chemical resistance. Autoclavable.



Clear

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
10175184	-	BWP0030PN	30	22	69	32	10
10136494	-	BWP0060PN	60	22	81	40	10
10728795	-	BLP0125P	125	32	104	50	10
10575105	-	BLP0150P	150	32	116	50	10
10287764	-	BWP0250P	250	37	140	60	10
10635325	-	BWP0500P	500	46	161	75	10
10034604	-	BWP1000P	1,000	55	205	90	5
10799145	-	BWP2000P	2,000	55	242	120	5

Amber

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
12094687	-	BWP0060AP	60	22	81	40	10
10739185	-	BWP0150AP	150	32	116	50	10
10329934	-	BWP0250AP	250	37	140	60	10
10554925	-	BWP0500AP	500	46	161	75	10
10799515	-	BWP1000AP	1,000	55	205	90	5
10267824	-	BWP2000AP	2,000	55	242	120	5

Have you thought about.....

Scoops, polypropylene, Fisherbrand®



Catalogue No	Legacy No	Capacity, mL	Length, overall, mm	Pack qty
11567852	-	25	135	5
11577852	-	50	160	5
11587852	-	100	200	5
11597852	-	250	260	1
11507862	-	500	315	1
11517862	-	1,000	385	1



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.