

Bottles, HDPE, wash, spray, adjustable, Fisherbrand®

Fisherbrand®

For hot, cold, aqueous, and solvent based solutions.

- Set to fine mist for TLC plates, position it to needle to rinse glassware, and pump a continuous, forceful jet on the big stream position for washing filters
- Instant off trigger release eliminates messy drips
- Chemical resistant
- Can be gas or chemically sterilised and cleaned in the laboratory dishwasher
- Zero contamination

The double headed 1,000mL sized bottle has a side fill opening to help maintain purity of solutions.

Catalogue No	Legacy No	Description	Pack qty
11795873	-	Double headed adjustable spray wash bottle, 1,000mL	1
11714356	-	Adjustable spray wash bottle, 240mL	3



Bottles, polyethylene, wash, trigger operated

With adjustable nozzle which produces a powerful intermittent or continuous jet or fine spray, head is mounted on a polyethylene bottle shaped to fit the hand.

Delivers 1.4mL per stroke.

Catalogue No	Legacy No	Capacity, mL
12413070	-	500



Bottles, HDPE, wash, venting, solvent labelled, Easy-Squeeze, Fisherbrand®



Fisherbrand®

For dispensing water and common solvents.

- Labelled and colour coded for most commonly used solvents
- Translucent with coloured closure
- High density polyethylene
- Comes with vented 38mm closure to prevent solvent drips
- 500mL capacity

Catalogue No	Legacy No	Label	Closure colour	Pack qty
11745233	-	Acetone	Red	6
11755233	-	Bleach	White	6
11765233	-	Deionised water	Natural	6
11775233	-	Distilled water	Natural	6
11785233	-	Ethanol	Natural	6
11795233	-	Isopropanol	Yellow	6
11705243	-	Methanol	Green	6
11715243	-	Saline	Natural	6
11725243	-	Soap	Natural	6
11735243	-	Universal	Natural	6
11745243	-	Water	Natural	6
11755243	-	Solvent	Natural	6
11765243	-	Cleaning	Natural	6



All bottles labelled for intended usage with specific hazardous substances are also labelled with the appropriate GHS hazard warning pictograms.