



Burettes, digital, Top Buret

eppendorf

- Continuous pulse-free dispensing technique for rapid titration
- Robust digital burette with fatigue free display and control elements
- TopBuret M dispensing rate 2,500µL per rotation
- TopBuret H dispensing rate 5,000µL per rotation
- Patented safety valve for de-airing without loss of reagent
- Dispensing range from 0.01mL to 999.9mL
- Telescopic aspiration tube
- Excellent chemical resistance
- Supplied with safety valve, telescopic operation tube, adjustable discharge tube, three adapters for 4mm, 38mm and 32mm necks and two 1.5V batteries

Catalogue No	Legacy No	Alt. No	Description
10236721	-	4965000017	TopBuret M
10236671	-	4965000025	TopBuret H

Accessory

Catalogue No	Legacy No	Alt. No	Description
10755972	-	4960851000	Dry tube



Burettes, digital

stuart

- Digital burette designed to combine speed and ease of use with high accuracy
- Display rotates through 360° so that it can always be positioned facing the user with the dispensing tube safely pointing away
- Large clear LCD display
- Accuracy $\pm 0.2\%$, reproducibility $\pm 0.01\%$, subdivision 0.01mL
- Supplied with calibration certificate
- Zero reset for multiple titrations
- Long life battery with power for 60,000 titrations
- Automatic shut off without losing reading to conserve power
- Simple disassembly for easy cleaning with autoclavable liquid pathway
- With BioCote® antimicrobial protection
- 30mm thread supplied with adapters for bottles with 33mm, 38mm and 45mm necks
- Supplied complete with works certificate but without reservoir

Catalogue No	Legacy No	Alt. No	Capacity, mL
11850501	-	BT30	30
11860501	-	BT50	50

Reservoir base with side neck, amber borosilicate

Catalogue No	Legacy No	Alt. No	Capacity, mL
11890501	-	DB1400A	1400
12812134	-	DB2500A	2500

Adapters

Catalogue No	Legacy No	Alt. No	Description
12822134	-	BTA40	For 40mm thread
12832134	-	BTA33	For 33mm thread
12842134	-	BTA38	For 38mm thread
12852134	-	BTA45	For 45mm thread