

Centrifuge tubes

Microcentrifuge tube cutter



This practical tool safely and easily removes the tips from plastic microcentrifuge tubes without the hazards of using razor blades or knives. Radioactive pellets caught in the tips may be retrieved without re-suspending them or risking splashed radioactive substances.



The cutter works well with tubes up to 1.5mL and composed of polyethylene, polypropylene or polyallomer. Not recommended for polystyrene. Maximum cut diameter is 9.4mm.

Catalogue No	Legacy No	Description
11730635	-	Microcentrifuge tube cutter
12395759	-	Replacement cutting blade

Tube racks, microcentrifuge, polypropylene, 4-way, Fisherbrand®



Four way heavy duty racks feature a unique system of tabs and slots which when linked together create a large rack.



Each rack can hold 4 x 50mL conical tubes, 12 x 15mL conical tubes, 32 x 1.5mL microtubes or 32 x 0.5mL microtubes. Rack measures 174mm x 95mm x 52mm. Autoclavable.

Catalogue No	Legacy No	Description	Pack qty
11750334	-	Microtube rack, 4-way, natural	5
11770045	-	Microtube rack, 4-way, blue	5
11760334	-	Microtube rack, 4-way, green	5
11780045	-	Microtube rack, 4-way, pink	5
11770334	-	Microtube rack, 4-way, yellow	5
11790045	-	Microtube rack, 4-way, orange	5
11700055	-	Microtube rack, 4-way, assorted colours (blue, green, pink, yellow and orange)	5

Tube rack, microcentrifuge, polypropylene, 3-way, Fisherbrand®



Three way polypropylene racks connect to one another on either end. Convenient finger slots are moulded onto the rack and each of the three tiers accommodates a different size tube. Lowest level has 24 wells to hold individual 0.2mL tubes or strips of 8 or 12 tubes. The middle tier has 14 wells to hold 0.5mL tubes. The top tier has 12 wells to hold 1.5mL tubes. Autoclavable. Rack measures 164mm x 123mm x 62mm and is supplied in assorted colours (blue, green, purple, yellow, and orange).




Catalogue No	Legacy No	Description	Pack qty
11928084	-	Microtube rack 3-way, assorted colours	6

Tube racks, microcentrifuge, polypropylene, 20 well, Fisherbrand®



Polypropylene racks for storing 20 tubes. Racks lock together on each side.



Tubes are easily removed using 'hold-it' finger slot on the front of the rack. Racks measure 213mm x 90mm x 50mm. Autoclavable.

Catalogue No	Legacy No	Description	Pack qty
11784016	-	Microtube rack 20 well, blue	6
11760045	-	Microtube rack 20 well, red	6
11704026	-	Microtube rack 20 well, white	6

Tube racks, microcentrifuge, polypropylene, cube, Fisherbrand®



User friendly cube racks.

Slide the polypropylene racks together in multiple orientations - work with different sized tubes at the same time. Each rack can hold 4 x 50mL tubes, 10 x 15mL conical tubes, 12 x 75mL or 100mL tubes or 16 x 1.5mL to 2.0mL tubes. Autoclavable. Racks measure 108mm x 113mm x 113mm.



Catalogue No	Legacy No	Description	Pack qty
11784726	-	Cube rack, assorted colours (one each of blue, green, yellow, orange, and pink)	5
11794726	-	Cube rack, natural	5

Tube racks, microcentrifuge, polypropylene, 80 well, Fisherbrand®



Polypropylene 80 well microtube racks ideal for freezer storage.

Racks hold 1.5mL to 2mL microtubes in a 5 x 16 array. Tubes sit in large wells with conical bottoms. Racks measure 225mm x 67mm x 28mm. Autoclavable.



Catalogue No	Legacy No	Description	Pack qty
11720344	-	Natural	5
11968074	-	Blue	5
11978074	-	Green	5
11988074	-	Pink	5
11998074	-	Yellow	5
11908084	-	Orange	5
11918084	-	Red	5
11728084	-	Assorted colours (one of each colour)	5