

Syringes, microvolume, fixed needle, 5µL to 500µL



Bevel tip needles listed, other lengths, gauges and tip styles available. Plunger protection syringes have an extended protection tube from the top of the flange to prevent plunger bending. Guided plunger syringe has extended barrel to guide plunger during injection and prevent bending. Repeating adapter is used to preset the syringe plunger to desired volume for reproducible dispensing.

Plunger protection

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10518502	-	001000	5	26	50
10607171	-	002000	10	26	50
10215322	-	002003	10	26	70

Catalogue No	Legacy No	Alt. No	Description
10127070	-	002030	Fixed needle, 6 pack of 10607171
12937536	-	002033	Fixed needle, 10 pack of 10607171

Guided plunger

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10442171	-	001400	5	26	50
10472171	-	002400	10	26	50

Standard plunger

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
11941134	-	003000	25	25	50
10229010	-	004000	50	25	50
10625471	-	005000	100	25	50
10559072	-	006000	250	25	50
10771422	-	007000	500	25	50

Syringes, microvolume, removable needle, 5µL to 500µL

- Bevel tip needles supplied



Other lengths, gauges, and tip styles available, please contact Customer Service.

Plunger protection syringes

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10025590	-	001050	5	26	50
10278720	-	002050	10	26	50

Catalogue No	Legacy No	Alt. No	Description
10288720	-	002080	6 pack of 10278720

Guided plunger

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10655851	-	001450	5	26	50
10452171	-	002450	10	26	50

With repeating adapter

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10025600	-	001090	5	26	50
10402001	-	003090	25	25	50
10177020	-	004090	50	25	50
10035470	-	005090	100	25	50
10472001	-	007090	500	25	50
10657171	-	002090	10	26	50
12997526	-	006090	250	25	50

Guided plunger and repeating adapter

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10646041	-	001490	5	26	50
10616231	-	002490	10	26	50

Standard plunger

Catalogue No	Legacy No	Alt. No	Volume, µL	Needle gauge	Needle length, mm
10065620	-	003050	25	25	50
10750852	-	004050	50	25	50
10228820	-	005050	100	25	50
10645851	-	006050	250	25	50
10421031	-	007050	500	25	50

