

Injection port liners, GC consumables, Thermo Scientific



Thermo
SCIENTIFIC

- Liners for Thermo Scientific and Agilent instruments
- Highly deactivated
- High reproducibility



Liners for TRACE™ and FOCUS™ SSL instruments

Catalogue No	Legacy No	Alt. No	Description	I.D. x O.D., mm	Length, mm	Pack qty
10651324	-	45350030	Split straight liner	5 x 8	105	5
10572075	-	45350031	Split straight liner	5 x 8	105	5
10197594	-	453T2131	Split straight liner with CarboFrit	5 x 8	105	5
10190553	-	45350032	Splitless straight liner	5 x 8	105	5
10455173	-	45350033	Splitless straight liner	5 x 8	105	5
10390855	-	453T2130	Splitless straight liner with CarboFrit	5 x 8	105	5
10416113	-	453T1905	Split FocusLiner for 50mm needle with quartz wool	5 x 8	105	5
10314273	-	453T2895	Split FocusLiner for 70mm needle with quartz wool	5 x 8	105	5
10751957	-	453T2999	Splitless FocusLiner for 50mm needle with quartz wool	5 x 8	105	5

Liners for Agilent instruments

Catalogue No	Legacy No	Alt. No	Description	I.D. x O.D., mm	Length, mm	Pack qty
10415923	-	453A2265	Split straight liner with quartz wool	4.0 x 6.3	78.5	5
10533615	-	453A1295	Split straight liner	4.0 x 6.3	78.5	5
10253753	-	453A1255	Split/splitless FocusLiner with quartz wool	4.0 x 6.3	78.5	5
10545435	-	453A1315	Split/splitless FocusLiner with single taper	4.0 x 6.3	78.5	5
10140893	-	453A1345	Split/splitless liner with single taper	4.0 x 6.3	78.5	5
10555435	-	453A1355	Split/splitless liner with double taper	4.0 x 6.3	78.5	5
10418484	-	453A1305	Split/splitless liner with recessed gooseneck, quartz wool	4.0 x 6.3	78.5	5
10055694	-	453A2375	Split/splitless FAST FocusLiner with single taper, quartz wool	2.3 x 6.3	78.5	5

Technical Support

For further information on any of the featured products contact our Technical Support Team



+44 1509 555939
www.fisher.co.uk
uk.international@thermofisher.com