

GC Accessories, Varian



Injectors: 1093 and 1094

Catalogue No	Legacy No	Alt. No	Type	I.D., mm	O.D., mm	Length, mm	Pack qty
10635841	-	092027	Restriction, SPI 0.25mm	0.5	4.6	54	5
10107210	-	092030	Restriction, SPI 0.25mm	0.8	4.6	54	5
10353748	-	092034	Restriction, SPI 0.5mm for 0.53mm I.D. on-columns	0.8	4.6	54	5
11942637	-	0726532	Viton O-ring for 1177	-	-	-	10
10113022	-	0726217	Graphite sealing ring for 1078 and 1079 injectors	-	-	-	10

Ferrules

15% Graphite/85% vespel


Injectors and detectors at atmospheric pressures

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Ferrule I.D., mm	Pack qty
10340611	-	072663	0.10 to 0.25	0.4	10
10300661	-	072654	0.32	0.5	10
10507922	-	072655	0.53	0.8	10

SilTite™ starter kits with 10 ferrules and 2 nuts

GC-MS Interface

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Ferrule I.D., mm
12634945	-	073300	0.10 to 0.25	0.4
12937101	-	073301	0.32	0.5
12947101	-	073302	0.53	0.8



Cutter, for metal tubing

VICI Jour®

The ideal tool for cutting 1.59mm (1/16 inch) stainless steel tubing and GLT. It makes smooth, uniform cuts that are suitable for zero dead volume connections without the need for de-burring or reaming.

Catalogue No	Legacy No	Alt. No	Description
12772680	-	JR-792	HPLC Tube cutter
12782680	-	JR-793	Replacement wheel

