

Containers, polypropylene, snap cap, 200mL to 960mL



- Available in three sizes: 200mL, 400mL and 960mL with or without a polyethylene leak proof snap cap
- 200mL has graduations at 25mL, 400mL at 50mL, 960mL is not graduated
- Choice of natural or white polypropylene in a variety of styles
- Sterile containers are indicated and are sterile by irradiation with peel away certificate of sterilisation on the outer carton

N.B.11733174 is not CE+IVD

200mL

Catalogue No	Legacy No	Description	Sterile, Yes/No	Colour	Pack qty
11819002	-	Container only	N	Natural	660
11871091	-	Container with separate cap	N	Natural	660
11730084	-	Cap on container	N	Natural	220
11838992	-	Cap on container	Y	Natural	220
11809012	-	Container with separate cap	N	White	660
11765085	-	Cap on container	Y	White	220

400mL

Catalogue No	Legacy No	Description	Sterile, Yes/No	Colour	Pack qty
11703174	-	Container only	N	Natural	460
11783164	-	Container with separate cap	N	Natural	460
11713174	-	Cap on container	N	Natural	185
11793164	-	Cap on container	Y	Natural	185
11743174	-	Container with separate cap	N	White	460
11753174	-	Cap on container	Y	White	185

960mL

Catalogue No	Legacy No	Description	Sterile, Yes/No	Colour	Pack qty
11783154	-	Container with separate cap	N	Natural	250
12612225	-	Cap on container	N	Natural	120
11703164	-	Cap on container	Y	Natural	80
11723174	-	Container with separate cap	N	White	250
12612205	-	Cap on container	N	White	120
11733174	-	Cap on container	Y	White	80



Containers, polypropylene, sterile, aseptic production, hinged cap, Fisherbrand®

Fisherbrand®

Suitable for sampling, carriage, and preservation of liquid and solid products.

The hinged cap, designed for one handed operation and to avoid drips, removes the risk of contamination from a cap misplaced during sampling.

Straight sided, disposable containers. Frosted area for writing.

11793144 and 12622135 taper slightly top to bottom (top diameter 75mm, bottom diameter 65mm) and have graduations in 25mL increments to 275mL.
'Biologically clean'.

Catalogue No	Legacy No	Capacity, mL	Diameter, mm	Height, mm	Colour	Pack qty
11713154	-	45	30	85	Clear	650
12672125	-	50	33	80	Blue	650
11743154	-	90	45	80	Natural	350
11753154	-	90	45	80	Blue	350
11793144	-	300	75/65	92	Natural	240
12622135	-	300	75/65	92	Blue	240

