

Cylinders, PMP, transparent, spouted, moulded graduations



Azlon

- Calibrated to ISO 6706, Class B tolerances
- Steam sterilisation at 121°C does not affect accuracy
- Easy to read moulded graduations, in transparent PMP
- Stable hexagonal bases with non-slip anti-suctions feature
- S.P.I. recycling codes and batch code on each cylinder

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
11323584	-	CT0010P	10	0.2	10
11333584	-	CT0025P	25	0.5	5
10566455	-	CT0050P	50	1.0	5
10665705	-	CT0100P	100	1.0	5
12070003	-	CT0250P	250	2.0	1
10218024	-	CT0500P	500	5.0	1
10524355	-	CT1000P	1,000	10.0	1
11363584	-	CT2000P	2,000	20.0	1



Cylinders, PPCO, graduated, spouted, Fisherbrand®



Fisherbrand

- Large rounded bases ribbed for reinforcement and stability
- Calibrated "To Contain/To Deliver" at 20°C - meets ASTM laboratory standards
- Chemical, heat and impact resistant makes these cylinders excellent for long term use
- Generous pour spout allows for easy transfer from cylinder to alternate containers
- No meniscus to confuse readings - eliminates guesswork
- Large base prevents rolling and tipping

PPCO cylinders

Catalogue No	Legacy No	Capacity, mL	Pack qty
11947884	-	10	12
11957884	-	25	12
11967884	-	50	12
11997874	-	100	8
11907884	-	250	8
11967874	-	500	6
11937874	-	1,000	4
11927874	-	2,000	2



Have you thought about.....

Flasks, borosilicate glass, volumetric, rimmed neck, Class A, ISO 1402, DIN 12664, Fisherbrand®



With one mark graduation, rim, polyethylene snap cap and blue printing.

Catalogue No	Legacy No	Capacity, mL	Neck I.D., mm	Neck size	Tolerance, mL	Pack qty
12902330	-	5	7	7/16	0.025	5
11577842	-	10	7	7/16	0.025	5
12922330	-	20	9	10/19	0.040	5
12932330	-	25	9	10/19	0.040	5
11507852	-	50	11	12/21	0.060	5
12952330	-	100	13	12/21	0.100	5



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.