

Rotary evaporator accessories, Rotavapor®



6

Evaporating flasks, pear shaped. With 29/32 joint

Catalogue No	Legacy No	Alt. No	Description
12732335	-	000431	Evaporating flask, pear shaped, 50mL
12582276	-	000432	Evaporating flask, pear shaped, 100mL
12732485	-	000433	Evaporating flask, pear shaped, 250mL
12381559	-	000434	Evaporating flask, pear shaped, 500mL
12782445	-	000435	Evaporating flask, pear shaped, 1,000mL
12792445	-	000436	Evaporating flask, pear shaped, 2,000mL
12702455	-	000437	Evaporating flask, pear shaped, 3,000mL

Powder drying flask, pear shaped with indentations to assist the drying of powdery substances. With 29/32 joint

7

Catalogue No	Legacy No	Alt. No	Description
12742335	-	000452	Powder flask, pear shaped with indentations, 500mL
12752335	-	000453	Powder flask, pear shaped with indentations, 1,000mL
12762335	-	000454	Powder flask, pear shaped with indentations, 2,000mL

Beaker flasks gross volume 500mL, net volume 150mL at 30° inclination. With 29/32 joint

Catalogue No	Legacy No	Alt. No	Description
12772335	-	034764	Beaker flask for liquids
12792335	-	034767	Beaker flask for powders, with indentations

Beaker flasks gross volume 1,500mL, net volume 450mL at 30° inclination. With 29/32 joint

Catalogue No	Legacy No	Alt. No	Description
12712345	-	034230	Beaker flask for liquids
12732345	-	034269	Beaker flask for powders, with indentations

Adapters

Catalogue No	Legacy No	Alt. No	Description
10128981	-	36576	Reitmeyer adapter

Receiving flasks. With spherical 35/20 ball joint

Catalogue No	Legacy No	Alt. No	Description
12712335	-	000422	Receiving flask with spherical 35/20 ball joint, not P+G coated, 100mL
12722335	-	000423	Receiving flask with spherical 35/20 ball joint, not P+G coated, 250mL
12772455	-	025264	Receiving flask with spherical 35/20 ball joint, P+G coated, 500mL
10609665	-	020728	Receiving flask with spherical 35/20 ball joint, P+G coated, 1,000mL
10505404	-	025625	Receiving flask with spherical 35/20 ball joint, P+G coated, 2,000mL
12782455	-	025266	Receiving flask with spherical 35/20 ball joint, P+G coated, 3,000mL

Receiving flasks, Low Temperature (LT) for C and CR cold trap assemblies. With spherical 35/20 ball joint

Catalogue No	Legacy No	Alt. No	Description
10601891	-	40774	Receiving flask, LT, P+G coated, 500mL
12792455	-	40775	Receiving flask, LT, P+G coated, 1,000mL
12786174	-	40776	Receiving flask, LT, P+G coated, 2,000mL
12712455	-	003275	Ball joint clamp for Büchi receiving flask

Flask support

Catalogue No	Legacy No	Alt. No	Description
12787173	-	48618	Flask support ring