



Filter holders, membrane holder, glass, 47mm, vacuum



For vacuum filtration of aqueous, organic or corrosive liquids for particulate contamination analysis.

Also recommended for HPLC solvent filtration. Vacuum connection to integral base and flask cap is above filtrate exit level, avoiding risk of drawing filtrate droplets into vacuum tubing. Laminated filters will not seal correctly in this holder.

Warning: Not recommended for use with flammable liquids.

Catalogue No	Legacy No	Alt. No	Description
10658451	-	XX1504700	With funnel, base, cap, clamp and fritted glass filter support

Spare parts

Catalogue No	Legacy No	Alt. No	Description
10410781	-	XX1004704	Funnel, 300mL, ground glass seal
10389730	-	XX1504702	Ground glass base and tubulated cap
10094440	-	XX1504705	Vacuum filtering flask, 1L
10596702	-	XX2004708	Stainless steel support screen, 47mm



Filter holders, membrane holder, glass, 90mm, vacuum



For applications where a 47mm diameter filter does not provide enough surface area, and an increase in speed of the filtering process or throughput is required.

The new 90mm glass filter holder provides vacuum filtering of liquids for analysis of particulate or biological contamination retained on the filter surface, or to produce medium quantities of ultra-clean filtrate. Sterile procedures require aseptically installing or sterile filter in the separately sterilised filter hold. Requires stainless smooth-tip filter forceps, a ground glass joint filtering flask and a vacuum source.

Warning: Not recommended for use with flammable liquids.

Catalogue No	Legacy No	Alt. No	Description
10085441	-	XX1009020	Glass filter holder with stainless steel screen, 90mm

Accessories

Catalogue No	Legacy No	Alt. No	Description
10094440	-	XX1504705	Vacuum filtering flask, 1L
10545292	-	XX1604706	Vacuum filtering flask, 5L, conical bottom
10419850	-	XX6602500	Solvent filtering dispenser
10217870	-	XX7100004	$\frac{3}{16}$ inch silicone tubing for vacuum use

Spare parts

Catalogue No	Legacy No	Alt. No	Description
10779201	-	XX1009003	Spring clamp, 90mm, anodised
10345141	-	XX1009000	Funnel, 90mm, glass, 1L
10015401	-	XX1009004	Glass tabulated base, 90mm
10769201	-	XX1009002	Support screen, PTFE coated, 90mm
10791073	-	XX1009010	PTFE gasket, pack of 10