

Filtration

filter units, vacuum



Filter units, bottle top, vacuum, with receiver



CORNING

Corning vacuum filter systems are designed for vacuum filtration of cell culture solutions or other aqueous solutions.

- New square membrane funnel design provides improved flow rates
- Stable profile receiver bottle with easy-grip sides and angled hose adapter minimises tipping
- Receiver converts to a storage bottle by simply unthreading filter and applying the sterile plug seal cap
- Filter funnel, dust cover and receiver/storage bottles are polystyrene
- Tubing adapter, filter adapter ring and plug seal cap are polyethylene
- Pore size and membrane type are printed on each unit
- Gamma radiation sterilised; certified non-pyrogenic

Choose from four membrane types: polyethersulfone (PES), cellulose acetate (CA), cellulose nitrate (CN) or nylon (NY). PES filters have the lowest protein binding and fastest flow rates.

Filter units - filter funnel, filter receiver and leak proof closure

Catalogue No	Legacy No	Alt. No	Membranes	Colour	Funnel, mL	Receiver capacity, mL	Pore size, μm	Membrane surface area, cm^2	Pack qty
10504211	-	431153	PES	Yellow	150	150	0.22	13.6	12
10609141	-	431096	PES	Yellow	250	250	0.22	19.6	12
10268951	-	431097	PES	Yellow	500	500	0.22	33.2	12
10296191	-	431098	PES	Yellow	1,000	1,000	0.22	54.5	12
10005341	-	431154	CA	Orange	150	150	0.22	13.6	12
10468181	-	430767	CA	Orange	250	250	0.22	19.6	12
10769922	-	430769	CA	Orange	500	500	0.22	33.2	12
10689311	-	430517	CA	Orange	1,000	1,000	0.22	54.5	12
10174191	-	431155	CA	Orange	150	150	0.45	13.6	12
10350412	-	430768	CA	Orange	250	250	0.45	19.6	12
10229051	-	430770	CA	Orange	500	500	0.45	33.2	12
10730693	-	430516	CA	Orange	1,000	1,000	0.45	54.5	12
10016110	-	430758	CN	Orange	500	500	0.22	33.2	12
10164191	-	430186	CN	Orange	1,000	1,000	0.22	54.5	12
10377941	-	430771	NY	Red	250	250	0.2	19.6	12
10106961	-	430773	NY	Red	500	500	0.2	33.2	12
10660502	-	430515	NY	Red	1,000	1,000	0.2	54.5	12

Choose online ordering

- Browse all our products and services with only a few clicks
- Take advantage of promotions throughout the year
- Orders are instantly and automatically processed
- Keep track of your budget using the order history functionality
- Access your contract prices
- Track your order online

UK and International: www.fisher.co.uk | **Ireland:** www.ie.fishersci.com

Sweden: www.fishersci.se | **Norway:** www.fishersci.no

Finland: www.fishersci.fi | **Denmark:** www.fishersci.dk

