



Filter units, flask top, vacuum, with receiver, Stericup™



- Tapered flask design for easy handling
- Recessed receiver bottom for easy stacking and storage of capped flasks
- Filter diameter: 73mm
- Funnel tab secures accessory AP20 glass microfibre prefilter (10258330 sold separately, please see our website) for hard to filter solutions
- Each unit is gamma radiation sterilised

Millipore® Express™ PES

For fast flow, low protein binding sterilisation of aqueous solutions.

Durapore® PVDF

For lowest protein binding sterilisation/clarification of biological solutions.

10411741, 10176660, 12457896, 12467896, 10440011: These products have been tested and validated for use in stem cell research applications. To determine their effects on mouse stem cell growth or differentiation, three lots of Stericup™-GP devices were used to filter media with LIF. Once filtered, this media was used to passage mouse stem cells five times to verify that Stericup™-GP filtration does not impact pluripotency of mouse stem cells.

Filter units - 150mL flask top filter, 150mL flask and cap

Catalogue No	Legacy No	Alt. No	Pore size, µm	Membranes	Pack qty
10411741	-	SCGPU01RE	0.22	PES	12
10626921	-	SCGVU01RE	0.22	PVDF	12
10722731	-	SCHVU01RE	0.45	PVDF	12

Filter units - 250mL flask top filter, 250mL flask and cap

Catalogue No	Legacy No	Alt. No	Pore size, µm	Membranes	Pack qty
10182431	-	SCVPU02RE	0.1	PES	12
10176660	-	SCGPU02RE	0.22	PES	12
10186660	-	SCGVU02RE	0.22	PVDF	12
10742731	-	SCHVU02RE	0.45	PVDF	12

Filter units - 500mL flask top filter, 500mL flask and cap

Catalogue No	Legacy No	Alt. No	Pore size, µm	Membranes	Pack qty
12457896	-	SCGPU05RE	0.22	PES	12
10156040	-	SCGVU05RE	0.22	PVDF	12
10412701	-	SCHVU05RE	0.45	PVDF	12

Filter units - 500mL flask top filter, 1,000mL flask and cap

Catalogue No	Legacy No	Alt. No	Pore size, µm	Membranes	Pack qty
12467896	-	SCGPU10RE	0.22	PES	12
10074360	-	SCGVU10RE	0.22	PVDF	12

Filter units - 1,000mL flask top filter, 1,000mL flask and cap

Catalogue No	Legacy No	Alt. No	Pore size, µm	Membranes	Pack qty
10440011	-	SCGPU11RE	0.22	PES	12
10196180	-	SCGVU11RE	0.22	PVDF	12
10196230	-	SCHVU11RE	0.45	PVDF	12

Replacement receiver flasks for Stericup™ filter units, with screw top lids

Catalogue No	Legacy No	Alt. No	Capacity, mL	Pack qty
10156280	-	SC00B02RE	250	12
10044410	-	SC00B05RE	500	12
10518822	-	SC00B10RE	1,000	12