

Flasks, Erlenmeyer, polycarbonate, sterile



CORNING

- Wide, easy access mouth
- Use as culture flasks or storage vessels
- Plug seal caps offer two positions, open for gas exchange or closed for liquid tight seal
- Vent cap supplied for use in open systems
- Plain bottom

Standard flasks

Catalogue No	Legacy No	Alt. No	Volume, mL	Graduations, mL	Neck size, mm	Cap style	Pack qty
10045411	-	431143	125	25	26	Vent	50
10022611	-	431144	250	25	31	Vent	50
10352742	-	431145	500	50	43	Vent	25
10114101	-	431147	1000	50	43	Vent	25
10654541	-	430421	125	25	26	Plug	50
10716851	-	430183	250	25	31	Plug	50
10664541	-	430422	500	50	43	Plug	25
10537041	-	431146	1000	50	43	Plug	25



Flasks, Erlenmeyer and Fernbach, polycarbonate, sterile



CORNING

Supplied with a flat, vented cap.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Style	Pack qty
10410482	-	431255	2,000	200	Erlenmeyer	6
10092591	-	431256	2,000	200	Erlenmeyer baffled bottom	6
10216821	-	431252	3,000	None	Fernbach	4
10410292	-	431253	3,000	None	Fernbach baffled bottom	4



Spinner flasks, glass, baffled, angled sidearms, ProCulture®

CORNING

- Baffles enhance aeration and agitation
- Unique impeller design for optimal stirring
- Sidearm design allows access for 25mL and 50mL pipettes
- Borosilicate glass vessels
- Flask centre cap

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre neck, mm	Sidearm neck, mm
10210532	-	4500-125	125	70	32
10678832	-	4500-250	250	70	32
10361972	-	4500-500	500	100	45
10477681	-	4500-1L	1,000	100	45
10371972	-	4500-3L	3,000	100	45
10158431	-	4502-3L	3,000	120	45
10648832	-	4500-6L	6,000	100	45
10542112	-	4502-6L	6,000	120	45

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre neck, mm	Sidearm neck, mm
10643243	-	4504-3L	3,000	140	45
10689022	-	4504-6L	6,000	140	45
10138481	-	4504-8L	8,000	140	45
10750383	-	4504-15L	15,000	140	45
10572492	-	4504-36L	36,000	140	45

