

# Funnels

glass, conical



## Funnels, soda lime glass, conical, Fisherbrand®

**Fisherbrand®**

Catalogue No	Legacy No	Top diameter, mm	Stem length, mm	Pack qty
11502423	-	35	40	1
11512423	-	45	45	1
11522423	-	50	50	1
11532423	-	55	55	1
11542423	-	60	60	1
11552423	-	70	70	1
11562423	-	100	100	1



## Funnels, borosilicate glass, conical, Fisherbrand®

**Fisherbrand®**

Catalogue No	Legacy No	Top diameter, mm	Pack qty
11572423	-	55	10
11582423	-	75	10
11592423	-	100	10
11502433	-	150	1
11512433	-	200	1




## Funnels, glass, conical

**PYREX®**

With short stem.

Catalogue No	Legacy No	Alt. No	Top diameter, mm	Stem length, mm	Stem O.D., mm
11337934	-	2140/01M	35	35	8
12923611	-	2140/02M	45	45	8
11357934	-	2140/03M	70	75	8
11367934	-	2140/06M	80	82	8
11377934	-	2140/10M	105	106	10
12099907	-	2140/11M	125	122	12
11397934	-	2140/12M	150	154	16
11307944	-	2140/14M	200	193	24



## Funnels, Duran® glass, conical, DIN 12445, ISO 4798

**SCHOTT  
DURAN**

With short stem.

Catalogue No	Legacy No	Alt. No	Top diameter, mm	Stem length, mm	Stem O.D., mm	Filter paper diameter, mm	Pack qty
13197164	-	213512308	35	35	6	45 to 55	10
13174594	-	213512805	45	45	6	55 to 70	10
13147434	-	213513304	55	55	8	70 to 90	10
13158023	-	213513801	70	70	8	110 to 125	10
13154604	-	213514103	80	80	10	125 to 150	10
13107164	-	213514609	100	100	10	150 to 185	10
13148033	-	213515708	150	150	16	240 to 270	10
13184254	-	213516104	200	175	26	320 to 385	1
13177244	-	213516901	300	175	30	500	1