

Homogenisers

Continuation

Homogeniser, stand fixed, Polytron® PT10-35 GT



Dispersing tools

Catalogue No	Legacy No	Alt. No	Shaft length, mm	Shaft dia., mm	Tip [rotor] dia., mm	Tip [rotor] speed [max], m/ s	Volume, working, mL	Vessel dia., mm	Typical applications	Price
10223442	-	PT-DA 3005/2EC	78	5.5	3	5.4	0.1 to 5	6	Dispersal of smallest amounts of plant, soft animal or soft human tissues in liquid e.g. for RNA extraction, can be treated in small conical Eppendorf® tubes	
10535613	-	PT-DA 3007/2EC	101	7.8	5	6.5	0.3 to 10	8	Dispersal of smallest amounts of plant, animal or human tissues in liquid e.g. for RNA extraction, breaking up to cells, dissolving solids, can be treated in round bottom Eppendorf® tubes	
10365212	-	PT-DA3012/2EC	154	12	9	12	2 to 250	8	Dispersal of smallest amounts of plant, animal or human tissues in liquid e.g. for RNA extraction, breaking up to cells, dissolving solids, can be treated in round bottom Eppendorf® tubes	
10263582	-	PT-DA3012/2WEC	154	12	9	12	2 to 250	8	Suitable for fibrous and tough sample material such as meat, fish, fruits, vegetables, etc.	
10435922	-	PT-DA 3012/2MEC	154	12	9	14	2 to 250	15	Dispersing plant, human, animal tissues, resins, pigments, extractions, production of suspension and emulsions	
10141582	-	PT-DA 3020/2EC	193	20	15	27	10 to 2,500	25	Dispersing plant, human, animal tissues, resins, pigments, extractions, production of suspension and emulsions	
10520373	-	PT-DA 3020/2MEC	193	20	15	27	10 to 2,000	25	Dispersing plant, human, animal tissues, resins, pigments, extractions, production of suspension and emulsions	
10530373	-	PT-DA3020/2WEC	193	20	15	27	10 to 2,500	25	Suitable for fibrous and tough sample material such as meat, fish, fruits, vegetables, etc	
13150673	-	PT-DA 3030/2G	250	30	22	30	100 to 1,000	35	Dispersing plant, human, animal tissues, resins, pigments, extractions, production of suspension and emulsions	
10345212	-	PT-DA30302WEC	250	30	22	30	100 to 4,000	35	Suitable for fibrous and tough sample material such as meat, fish, fruits, vegetables, etc.	
10530563	-	PT-DA 3030/4EC	250	30	26	30	100 to 4,000	35	Dispersing plant, human, animal tissues, resins, pigments, extractions, production of finest suspensions and emulsions, gassing, extractions	

Choose online ordering

- Browse all our products and services with only a few clicks
- Take advantage of promotions throughout the year
- Orders are instantly and automatically processed
- Keep track of your budget using the order history functionality
- Access your contract prices
- Track your order online

UK and International: www.fisher.co.uk | **Ireland:** www.ie.fishersci.com

Sweden: www.fishersci.se | **Norway:** www.fishersci.no

Finland: www.fishersci.fi | **Denmark:** www.fishersci.dk

