

Incubators, INB, INE and INP Series

INE series

Catalogue No	Legacy No	Alt. No	Capacity, L
10758864	-	INE200	32
10162692	-	INE300	39
10172692	-	INE400	53
10345952	-	INE500	108
10040792	-	INE600	256
10477722	-	INE700	416
10355952	-	INE800	749

INP series

Catalogue No	Legacy No	Alt. No	Capacity, L
10182592	-	INP200	39
10466552	-	INP300	53
10788474	-	INP400	256
10798474	-	INP500	416
10020752	-	INP600	749
10122692	-	INP700	416
10030752	-	INP800	749

INB series

Catalogue No	Legacy No	Alt. No	Capacity, L
10799064	-	INB200	32
10010982	-	INB300	39
10435982	-	INB400	53
10502363	-	INB500	108

Accessories

Catalogue No	Legacy No	Alt. No	Description
11510066	-	G7	Tube extension
11554136	-	V9	Additional Steri-Card (for sterilisation of working chamber) for INP

Shelves, perforated

Catalogue No	Legacy No	Alt. No	Description
11579975	-	E0-200	For 32L units
11589975	-	E0-300	For 39L units
11599975	-	E0-400	For 53L units
11509985	-	E0-500	For 108L units
11519985	-	E0-600	For 256L units
11529985	-	E0-700	For 416L units
11539985	-	E0-800	For 749L units

Trays, non-perforated

Catalogue No	Legacy No	Alt. No	Description
11559985	-	E2-200	For 32L units
11569985	-	E2-300	For 39L units
11579985	-	E2-400	For 53L units
11589985	-	E2-500	For 108L units
11599985	-	E2-600	For 256L units
11509995	-	E2-700	For 416L units
11519995	-	E2-800	For 749L units

Wall brackets

Catalogue No	Legacy No	Alt. No	Description
11590006	-	G0-200	For 32L units
11500016	-	G0-300	For 39L units
11510016	-	G0-400	For 53L units
11520016	-	G0-500	For 108L units
11530016	-	G0-600	For 256L units
11540016	-	G0-700	For 416L units

Sub frames

Catalogue No	Legacy No	Alt. No	Description
11550016	-	G5-500	For 108L units
11560016	-	G5-600	For 256L units
11570016	-	G5-700	For 416L units
11580016	-	G6-500	With castors, for 108L units
11590016	-	G6-600	With castors, for 256L units
11500026	-	G6-700	With castors, for 416L units

Factory Fitted Options - see side panel for ordering procedure - NOT AVAILABLE FOR INB SERIES

Catalogue No	Legacy No	Description
12795058	-	Adjustable temperature limiter, for example for heat-sensitive loads for INE series
10404533	-	OQ check list including temperature distribution survey to DIN 12 880 for nine measuring points with works test data for incubator as support for validation by customer
11921469	-	Entry port (left centre top)
10234238	-	Entry port (right centre top)
11931469	-	Standard Plus with water jacket
11941469	-	Right - other port (23mm dia.) in special position
11981429	-	Standard Plus with pre/post vacuum system (50mbarA)
11901439	-	Standard Plus (includes automatic water fill and condensate unit)
10571450	-	Flexible Pt100 for positioning in chamber or in load with socket, 4 pin, for external temperature recording, load temp
15134642	-	Potential free contact (24V/2A) with socket, 3 pin, for external temperature recording
11986578	-	Potential free contact for combined fault message (e.g., supply failure, sensor fault, fuse)
10217984	-	Floating triple contact, for signal generation, controlled by programme segment (using PC), for a total of three freely selected functions to be activated (e.g., activation of audible and visual signals, exhaust motors, fans, stirrers etc.)
11758907	-	IQ check list with works test data for humidity chamber as support for validation by customer
11921449	-	OQ check list including one free-selectable temperature distribution survey to DIN 12880 for 9 measuring points with works test data for humidity as support for validation by customer

Works calibration certificate at 37°C available, please enquire.