

## Lubricants

### WD 40 spray

Lubricant, 400mL, aerosol can.

### Silicone grease

Resistant to most chemicals. Stable over temperature range -40°C to 200°C. For high vacuum purposes. Must not be used with volumetric glassware because silicone contamination distorts the shape of the meniscus.

### PTFE, aerosol

A dry-film aerosol lubricant comprising a PTFE dispersion with a non-toxic bonding agent. For glass or metal joints. A useful anti-stick lubricant for glass joints. Thermally stable up to 275°C.

### Switch cleaner

Aerosol switch cleaner and lubricant.



| Catalogue No | Legacy No | Description     | Quantity    |
|--------------|-----------|-----------------|-------------|
| 12836824     | -         | WD 40 spray     | 415mL       |
| 12302058     | -         | Silicone grease | tube of 50g |
| 12856824     | -         | PTFE, aerosol   | 400mL       |
| 12865933     | -         | Switch cleaner  | 176mL       |

## Magnifiers, hand held, Fisherbrand®



**Fisherbrand®**

For precise scientific, plant and fieldwork.

- 11794386 and 11704396 illuminate focused images uniformly (supplied with AA batteries)
- 11714396 with adjustable focus



|   | Catalogue No | Legacy No | Description                       | Dimensions, mm       | Magnification            |
|---|--------------|-----------|-----------------------------------|----------------------|--------------------------|
| 1 | 11794386     | -         | Large illuminated magnifier       | 70 x 70 (window)     | 2x with 4x spot          |
| 2 | 11704396     | -         | Medium illuminated dual magnifier | 44.5 x 44.5 (window) | 3x with 10x spot         |
| 3 | 11714396     | -         | Small magnifier, not illuminated  | 25.44 (diameter)     | 8x with adjustable focus |

## Magnifier, illuminated, pocket microscope, portable, Fisherbrand®

**Fisherbrand®**

Illuminated pocket microscope ideal for inspection, quality control, field use, electronics, metal analysis, clinical observation, cloth/paper examination and plant studies. Adjustable light source directs the bright bulb on to the exact area to be examined. An adjustable focus control knob produces high quality imaging. Chemical and shock resistant ABS plastic case. Supplied with AA batteries and a snap on slide holder.

| Catalogue No | Legacy No | Description                               | Dimensions [w x d x h], mm |
|--------------|-----------|---|----------------------------|
| 11724396     | -         | Illuminated microscope, zoom, 60x to 100x | 140 x 38 x 23              |

### Spare part

| Catalogue No | Legacy No | Alt. No | Description                     | Pack qty |
|--------------|-----------|---------|---------------------------------|----------|
| 12747686     | -         | MN1500  | Batteries, alkaline, 1.5V (LR6) | 4        |



All products are for research use only, unless otherwise indicated.