



Microplates, 96 well, polystyrene, non-treated, Thermo Scientific Sterilin™



Thermo Scientific Sterilin™ clear multiwell plates are manufactured from premium grade virgin polystyrene and offer excellent optical clarity. Black plates are ideal for fluorescence work.

Suitable for use at 340nm. U wells hold 330µL, V wells hold 310µL and the flat based wells hold 400µL in volume. All Thermo Scientific Sterilin™ multiwell plates are now CE marked in accordance with the In Vitro Diagnostic Directive (98/79/EC).

All plates are non-treated. Non-sterile plates provide a medium binding surface. Irradiation process for the sterile plates makes the surface more hydrophilic, therefore generating a high binding product.

Catalogue No	Legacy No	Alt. No	Type	Colour	Sterile, Yes/No	Pack qty
11399153	-	612U96	U well	Clear	Yes	50
11309163	-	612V96	V well	Clear	Yes	50
11319163	-	612F96	Flat bottomed	Clear	Yes	50
11329163	-	642000	Lids	Clear	Yes	50
11339163	-	611U96	U well	Clear	No	50
11897403	-	611V96	V well	Clear	No	50
11349163	-	611F96	Flat bottomed	Clear	No	50
11359163	-	611F96BK	Flat bottomed	Solid black	No	50

All products are for research use only, unless otherwise indicated.



Dishes, polystyrene, sterile, CE marked for IVF, Thermo Scientific Nunc®



- Intended for IVF use
- Unique certification based on a 1-cell stage mouse embryo toxicity test (MEA test); the certificate confirms sterility (SAL 10⁻⁶), non-pyrogenicity and that the material has passed USP Class VI test
- The embryo toxicity test is a release test and the product will only be sold after it has passed the test
- Supplied with lids
- Manufactured under an ISO 13485 registered system to ensure consistent quality
- Ergonomic design for easier handling
- Flat bottomed for optimal heat transfer (ICSI = Intra-cytoplasmic sperm injection)

Catalogue No	Legacy No	Alt. No	Surface	No. of wells	Dimensions, mm	Culture area, cm²	Inner pack qty	Pack qty
10404532	-	144444	Nuncclon™Δ	4	66 x 66	1.9/well	4	120
10756864	-	179830	Non-treated	4	66 x 66	1.9/well	4	120
10170753	-	150255	Non-treated	1	40 dia. x 12	8.8	10	500
10262623	-	150270	Non-treated	1	60 dia. x 15	21.5	10	400
10021963	-	150360	Non-treated	1	90 dia. x 17	56.7	10	150
10713697	-	150260	Non-treated	1, centre well	16 x 55	-	1	120
10157524	-	150265	Non-treated	1, ICSI plate	9 x 51	-	1	120

All products are for research use only, unless otherwise indicated.

For our full range of microplates, please see our website

