

## Heating and drying ovens, Thermo Scientific Heratherm®


**Thermo**  
SCIENTIFIC

## Advanced Protocol Security

Catalogue No	Legacy No	Alt. No	Type	Volume, L	Dimensions, external [w x d x h], mm	Dimensions, internal [w x d x h], mm	Mass, kg	No. of shelves [max/standard]
10190010	-	51028145	Gravity convection, coated exterior	61	530 x 565 x 720	328 x 389 x 480	44	13/2
10100020	-	51028146	Gravity convection, coated exterior	99	640 x 565 x 820	438 x 389 x 580	55	16/2
10307539	-	51028147	Gravity convection, coated exterior	168	640 x 738 x 920	438 x 564 x 680	69	19/2
10222320	-	51028544	Gravity convection, stainless steel exterior	61	530 x 565 x 720	328 x 389 x 480	44	13/2
10232320	-	51028545	Gravity convection, stainless steel exterior	99	640 x 565 x 820	438 x 389 x 580	55	16/2
10110080	-	51028546	Gravity convection, stainless steel exterior	168	640 x 738 x 920	438 x 564 x 680	69	19/2
10100550	-	51028154	Mechanical convection, coated exterior	62	530 x 565 x 720	328 x 343 x 508	44	13/2
10820181	-	51028155	Mechanical convection, coated exterior	97	640 x 565 x 820	438 x 343 x 608	55	16/2
10570991	-	51028156	Mechanical convection, coated exterior	170	640 x 738 x 920	464 x 518 x 708	69	19/2
10660421	-	51028547	Mechanical convection, stainless steel exterior	62	530 x 565 x 720	328 x 343 x 508	44	13/2
10672341	-	51028548	Mechanical convection, stainless steel exterior	97	640 x 565 x 820	438 x 343 x 608	55	16/2
10775182	-	51028549	Mechanical convection, stainless steel exterior	170	640 x 738 x 920	464 x 518 x 708	69	19/2

## Accessories

Catalogue No	Legacy No	Alt. No	Description
13487669	-	50128881	Reinforced shelf for 180L models, max total capacity is 250kg
13497669	-	50128880	Reinforced shelf, 180L, gravity convection shelf with maximum capacity of 70kg for 180L models, gravity convection of Advanced Protocol and advanced protocol and advanced protocol security ovens; max total capacity is 250kg
10303990	-	50127761	Additional wire mesh shelf for 60L gravity convection models, including 2 shelf supports
10130250	-	50127762	Additional wire mesh shelf for 100L gravity convection models, including 2 shelf supports
10622131	-	50127763	Additional wire mesh shelf for 180L gravity convection models, including 2 shelf supports
10744412	-	50127764	Wire mesh shelf for IMH60, IMH60-S, IGS60, including 2 shelf supports
10754412	-	50127765	Wire mesh shelf for IMH100, IMH100-S, IGS100, including 2 shelf supports
10262150	-	50127766	Wire mesh shelf for IMH180, IMH180-S, IGS180, including 2 shelf supports
10784022	-	50126665	Stacking kit for 60L models - not required but recommended if decontamination cycle is performed in lower unit
10801311	-	50126666	Stacking kit for 100L models - not required but recommended if decontamination cycle is performed in lower unit
10019729	-	50126667	Stacking kit for 180L models - not required but recommended if decontamination cycle is performed in lower unit
10811311	-	50127741	Support stand with casters for 60L models
10304130	-	50127742	Support stand with casters for 100L models
10724602	-	50127743	Support stand with casters for 180L models
10869470	-	50127767	Sample sensor to measure and display exact sample temperature for Heratherm® Advanced Protocol Security models
12612785	-	50128887	Kit for under bench installation; for GP ovens; consists of exhaust tubing (1.5m length) and details on wall; maximum temperature for underbench installation is 250°C

## Factory Fitted Options; please order in conjunction with desired oven

Catalogue No	Legacy No	Description
11894413	-	Additional access port on centre of left side of unit; dia. 24mm for all incubators
11804423	-	Additional access port on centre of left side of unit; dia. 58mm for all incubators
11814423	-	Additional access port on centre of right side of unit; dia. 24mm for all incubators
11824426	-	Additional access port on centre of right side of unit; dia. 58mm for all incubators
11834423	-	Additional access port on centre of top of unit; dia. 24mm for all incubators
11844423	-	Additional access port on centre of top of unit; dia. 58mm for all incubators
13254779	-	Door hinge on left side. Available for General Protocol, Advanced Protocol and Advanced Protocol Security units. Coated exterior only


**Factory Fitted Options  
Ordering Procedure**

- 1) Specify catalogue number of option(s) required to be factory fitted
- 2) Place your order in the usual manner

The prices given include fitting. Should you have any queries please contact Customer Service.

**ATTENTION**  
**FACTORY FITTED OPTIONS MUST BE SPECIFIED AT THE TIME OF ORDERING THE MAIN UNIT**