

Absorbent bench paper, Fisherbrand®



This polyethylene backed 125g/m² absorbent paper is suitable to protect all types of laboratory work surfaces. Available in both sheet or roll form.

Catalogue No	Legacy No	Dimensions	Pack qty
11527362	-	570mm x 460mm sheets	100
11758722	-	460mm x 50m roll	1
11547362	-	920mm x 50m roll	1




Absorbent bench protector, Benchkote®



An absorbent, impermeable material designed to protect laboratory surfaces against hazardous spills. The material features a high quality smooth, absorbent Whatman paper which quickly absorbs liquid spills, and a laminated polyethylene layer that prevents flow through the work surface. After use, the sheet is incinerated or disposed of according to local regulations.

Catalogue No	Legacy No	Alt. No	Dimensions	Pack qty
10682941	-	2300-916	460mm x 570mm, sheets	50
11820815	-	2300-917	400mm x 570mm, sheets	100
11752187	-	2300-599	460mm x 570mm, sheets in pad	50
10672751	-	2300-731	460mm x 50m, roll	1
10008000	-	2300-772	920mm x 50m, roll	1



Benchkote® Plus thicker more absorbent material for more demanding applications and can absorb in excess of 0.75L

Catalogue No	Legacy No	Alt. No	Dimensions	Pack qty
10164021	-	2301-6150	500mm x 600mm	50
10144861	-	2301-6160	600mm x 50m, roll	1

Absorbent bench protector, BenchGuard™, Thermo Scientific Sterilin®



Plastic coated white paper, highly absorbent. Protects surfaces from accidental spillage of chemicals and biohazards.

Catalogue No	Legacy No	Alt. No	Absorbency	Dimensions
11359183	-	BG50	400mL/m ²	490mm x 50m roll
11379183	-	BG92	400mL/m ²	920mm x 50m roll

BenchGuard™ Extra

Catalogue No	Legacy No	Alt. No	Absorbency	Dimensions	Pack qty
11389183	-	BG50E	800mL/m ²	490mm x 50m roll	1
11399183	-	BG60E	800mL/m ²	490mm x 600mm sheets	50



Absorbent bench liners, Labmat™



Bench and drawer liner simplifies clean up chores and safeguards work surfaces from stains, spills and wear.

The flexible polyethylene film backing is non-slip and prevents seepage of most solutions. The highly absorbent upper surface traps water, dyes, corrosive chemicals and radioactive solutions. Each roll is supplied in a handy dispenser carton. Required length is easily cut off with scissors while the roll is retained in the dispenser.

Catalogue No	Legacy No	Alt. No	Description	Dimensions, mm
11367234	-	F24675-0000	Labmat™ roll	490mm x 50m roll
11342194	-	F24675-1000	Biohazard Labmat™	508mm x 15.25m

