

Pipettes, straight, serological, disposable glass, Pyrex®



CORNING

- Calibrated to deliver with blow out
- All pipettes have negative graduations
- Steriview wrappers (paper/plastic)
- Colours conform to ISO standards
- Long, slender, tapered tips to make pipetting go faster and easier

Individually wrapped, plugged, sterile

Catalogue No	Legacy No	Alt. No	Capacity, mL	Sleeve	Pack type	Pack qty
11303845	-	7077-5N	5	1	Individually wrapped	120
11313845	-	7077-10N	10	1	Individually wrapped	100
13053093	-	7078-10N	10	20	Bulk wrapped	120



Pipettes, straight, serological, polystyrene, Fisherbrand®



Fisherbrand®

Sterilised by electron beam treatment these high quality, blow out, plugged, crystal polystyrene, mono moulded, serological pipettes are particularly suitable for food microbiological applications. Marked with strong, black, printed graduations and negative graduations. Colour coded and available in free flow open ended style or conventional taper style, they are accurate to better than 1.5% tolerance.

Open ended style


Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Negative gradient, mL	Colour code	Inner pack qty	Pack qty
11542862	-	1	0.01	0.2	Yellow	25	1,000

Conventional taper, individually wrapped

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Negative gradient, mL	Colour code	Pack qty
11552862	-	1	0.01	0.2	Yellow	500
11562862	-	2	0.01	0.2	Green	500
11592862	-	25	0.2	5.0	Red	100

Conventional taper

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Negative gradient, mL	Colour code	Inner pack qty	Pack qty
11572862	-	2	0.01	0.2	Green	25	1,000
11582862	-	5	0.1	0.15	Blue	25	375



Pipettes, straight, serological, polystyrene, sterile, Thermo Scientific Nunc®



Thermo
SCIENTIFIC

nunc

- Accurate, disposable plastic pipettes plugged and sterilised
- Sterility assurance level (SAL) of 10⁻⁶
- Bold clear graduations with easy to read scale markings
- Colour coded packaging for ease in sorting and selecting the correct size
- Individually wrapped in paper peelable bags
- Convenient extra graduations to full pipette volume
- Non-pyrogenic

Catalogue No	Legacy No	Alt. No	Volume, mL	Subdivisions, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
10140771	-	159609	1	0.01	0.3	Black	200	1,000
10489472	-	159617	2	0.01	0.3	Black	125	500
10739731	-	159625	5	0.1	2	Blue	50	200
10292431	-	159633	10	0.1	3	Orange	50	200
10252531	-	159641	25	0.2	10	Green	50	200
10212581	-	159668	50	0.5	10	Purple	25	100

