



Pipettes, micro sampling



IntraMARK micropipettes, BLAUBRAND®, with ring mark conform to DIN ISO 7550 and are certified for volumes $\geq 5\mu\text{L}$, accuracy $\pm 0.25\%$, reproducibility $\geq 0.5\%$ (size $5\mu\text{L}$; accuracy $\pm 0.3\%$; reproducibility $\geq 0.6\%$). Fire polished surfaces inside and out are particularly smooth; the aspiration end is also fire polished to minimise pipette adapter wear. Calibrated to contain ring mark. Colour code to ISO for clear identification. Length 127mm.

Catalogue No	Legacy No	Alt. No	Capacity, μL	Colour	Pack qty
10450233	-	708707	1-5	White	250
10788181	-	708709	10	Orange	250
10788371	-	708718	20	Black	250
10354261	-	708722	25	2 x White	250
10252611	-	708733	50	Green	250
10653981	-	708744	100	Blue	250
12961868	-	708757	200	Red	250
11718208	-	708727	40	2 x Red	250
11728208	-	708728	20 + 40	2 x Red	250



Pipette bulbs, rubber, red



Contoured shape and heavy wall neck increase grip for easier use with pipettes and syringes.

The 150mL capacity bulb has reinforcing ribs for use with larger pipettes.

Rubber bulbs

Catalogue No	Legacy No	Capacity, mL	Neck diameter, mm
11967984	-	12	4
11977984	-	15	4
11987984	-	30	6
11997984	-	60	6
11314055	-	150	6



Pipette bulbs, rubber teat

Catalogue No	Legacy No	Bore size, mm	Capacity [nominal], mL	Pack qty
12406180	-	5	1	20
12416180	-	5	1.5	20
12426180	-	6	2.5	20
12436180	-	6	3.75	20

Catalogue No	Legacy No	Description	Pack qty
12446180	-	Mixed pack of 5 of each of the above	20