



Pipettor tips, Tecan® compatible



- TT-1000 tips made from conductive LLS black polypropylene
- ZT and ZTF tips made from clear polypropylene
- With or without aerosol resistant filters
- Liquid level sensing (TT-1000-CBK), standard or Maxymum Recovery®
- Non-sterile or sterile
- Certified DNase free, RNase free and pyrogen safe

Compatible with Tecan® Genesis Freedom®, Freedom Evo® and Miniprep with LIHA, liquid level sensing

	Catalogue No	Legacy No	Alt. No	Type	Volume, max., µL	For use with	Sterile, Yes/No	Inner pack qty	Pack qty
1	11364285	-	TT-1000-CBK-HTR	Conductive, non-filtered	1,000	LiHa-1, -4 and -8	No	96	1,536
	12521288	-	TT-1000-CBK-HTR-24	Conductive, non-filtered	1,000	LiHa-1, -4 and -8	No	96	4,608
	11384285	-	TT-1000-CBK-HTR-S	Conductive, non-filtered	1,000	LiHa-1, -4 and -8	Yes	96	1,536
	11354305	-	TTF-1000-CBK-HTR-S	Conductive, filtered	1,000	LiHa-1, -4 and -8	Yes	96	1,536

Compatible with Tecan® Freedom Evo® Temo and Aquarius

	Catalogue No	Legacy No	Alt. No	Type	Volume, max., µL	For use with	Sterile, Yes/No	Inner pack qty	Pack qty
2	12706689	-	ZT-50-R	Non-filtered	50	MCA96, TEM096	No	96	4,800
	12716689	-	ZT-50-R-S	Non-filtered	50	MCA96, TEM096	Yes	96	4,800
	12796769	-	ZTF-50-R-S	Filtered	50	MCA96, TEM096	Yes	96	4,800
	12786669	-	ZT-100-R	Non-filtered	100	MCA96, TEM096	No	96	4,800
	12796669	-	ZT-100-R-S	Non-filtered	100	MCA96, TEM096	Yes	96	4,800
	12736769	-	ZTF-100-R-S	Filtered	80	MCA96, TEM096	Yes	96	4,800
	12726679	-	ZT-200-R	Non-filtered	200	MCA96, TEM096	No	96	4,800
	12736679	-	ZT-200-R-S	Non-filtered	200	MCA96, TEM096	Yes	96	4,800
	12756769	-	ZTF-200-R-S	Filtered	180	MCA96, TEM096	Yes	96	4,800
	12756679	-	ZT-205-WB-R	Wide bore, non-filtered	200	MCA96, TEM096	No	96	4,800
	12766679	-	ZT-205-WB-R-S	Wide bore, non-filtered	200	MCA96, TEM096	Yes	96	4,800
	12776769	-	ZTF-205-WB-R-S	Wide bore, filtered	180	MCA96, TEM096	Yes	96	4,800

Compatible with Tecan® Freedom Evo® Temo and Aquarius, Maxymum Recovery®

	Catalogue No	Legacy No	Alt. No	Type	Volume, max., µL	For use with	Sterile, Yes/No	Inner pack qty	Pack qty
	11336064	-	ZT-50-L-R	Non-filtered	50	MCA96, TEM096	No	96	4,800
	12796679	-	ZT-50-L-R-S	Non-filtered	50	MCA96, TEM096	Yes	96	4,800
3	12786769	-	ZTF-50-L-R-S	Filtered	50	MCA96, TEM096	Yes	96	4,800
4	11346064	-	ZT-100-L-R	Non-filtered	100	MCA96, TEM096	No	96	4,800
	12776669	-	ZT-100-L-R-S	Non-filtered	100	MCA96, TEM096	Yes	96	4,800
5	12726769	-	ZTF-100-L-R-S	Filtered	80	MCA96, TEM096	Yes	96	4,800
6	11356064	-	ZT-200-L-R	Non-filtered	200	MCA96, TEM096	No	96	4,800
	12716679	-	ZT-200-L-R-S	Non-filtered	200	MCA96, TEM096	Yes	96	4,800
7	12746769	-	ZTF-200-L-R-S	Filtered	180	MCA96, TEM096	Yes	96	4,800
8	12746679	-	ZT-205-WB-L-R	Wide bore, non-filtered	200	MCA96, TEM096	No	96	4,800

