



Pipettor starter kits, Proline® Plus



Proline® Plus starter kits include pipettors (as listed below) and pipettor holder, calibration tool and pipettor tips for each pipettor and literature.

Proline® Plus kit 1

- Proline® Plus pipettor, 0.1µL to 3µL
- Proline® Plus pipettor, 0.5µL to 10µL

Proline® Plus kit 2

- Proline® Plus pipettor, 0.5µL to 10µL
- Proline® Plus pipettor, 10µL to 100µL
- Proline® Plus pipettor, 100µL to 1,000µL

Proline® Plus kit 3

- Proline® Plus pipettor, 2µL to 20µL
- Proline® Plus pipettor, 20µL to 200µL
- Proline® Plus pipettor, 100µL to 1,000µL

Proline® Plus kit 4

- Proline® Plus pipettor, 500µL to 5,000µL
- Proline® Plus pipettor, 1,000µL to 10,000µL

Catalogue No	Legacy No	Alt. No	Description
13114474	-	728650	Proline® Plus kit 1
13164394	-	728651	Proline® Plus kit 2
13194524	-	728652	Proline® Plus kit 3
13171604	-	728653	Proline® Plus kit 4



Pipettor starter kits, mLINE®



The Sartorius Biohit mLINE® starter kits offer an opportunity to test and get started with mLINE® pipettors, choose between two starter kits both containing three single channel mLINE® pipettors. Both kits also include colour coded Optifit tip racks for every pipettor and many useful accessories such as pipettor holders, colour coding caps and a calibration tool, which also acts as a tube opener.

mLINE® starter kit 1

- mLINE® pipettor 0.5µL to 10µL
- mLINE® pipettor 10µL to 100µL
- mLINE® pipettor 100µL to 1,000µL
- 3 x pipettor holders
- Colour coding caps
- Pen
- Optifit tips 0.2µL to 10µL, single tray
- Optifit tips 0.5µL to 300µL, single tray
- Optifit tips 10µL to 1,000µL, single tray
- Calibration tool/tube opener
- Literature

mLINE® starter kit 2

- mLINE® pipettor 2µL to 20µL
- mLINE® pipettor 20µL to 200µL
- mLINE® pipettor 100µL to 1,000µL
- 3 x pipettor holders
- Colour coding caps
- Pen
- 2 x Optifit tips 0.5µL to 300µL, single tray
- Optifit tips 10µL to 1,000µL, single tray
- Calibration tool/tube opener
- Literature

	Catalogue No	Legacy No	Alt. No	Description
1	13164664	-	725651	mLINE® starter kit 1
2	13167554	-	725652	mLINE® starter kit 2

Serving Science through our catalogue ranges:

- Chemicals
- Laboratory
- Life Science
- Safety



ONE SOURCE. INFINITE SOLUTIONS.®