


Continuation



Condensers



Jacketed coil, screwthread connection

| Catalogue No | Legacy No | Alt. No | Surface area, x10 ⁻² m ² | Socket size | Cone size | Length, effective, mm | Length, mm |
|--------------|-----------|-----------|--|-------------|-----------|-----------------------|------------|
| 1 12000073 | - | CX6/22/SC | 4.0 | 19/26 | 19/26 | 160 | 360 |
| 12010073 | - | CX6/33/SC | 4.0 | 24/29 | 24/29 | 160 | 365 |
| 12498836 | - | CX6/05/SC | 4.0 | - | 34/35 | 160 | 380 |

Liebig, ground cone and socket, screwthread connection

| Catalogue No | Legacy No | Alt. No | Surface area, x10 ⁻² m ² | Socket size | Cone size | Length, effective, mm | Length, mm |
|--------------|-----------|----------|--|-------------|-----------|-----------------------|------------|
| 2 12019652 | - | C1/11/SC | 0.6 | 14/23 | 14/23 | 160 | 260 |
| 11301855 | - | C1/12/SC | 1.0 | 19/26 | 19/26 | 200 | 340 |
| 12039652 | - | C1/13/SC | 1.0 | 24/29 | 24/29 | 200 | 390 |
| 12049652 | - | C1/22/SC | 2.0 | 19/26 | 19/26 | 400 | 530 |
| 11321855 | - | C1/23/SC | 2.0 | 24/29 | 24/29 | 400 | 535 |



Dean & Stark receivers, BS 756



Receivers for light entrainers

| Catalogue No | Legacy No | Alt. No | Capacity, mL | Graduations, mL | Socket | Cone size |
|--------------|-----------|-----------|--------------|-----------------|--------|-----------|
| 10422025 | - | WE6/23/2 | 2.0 | 0.05 | 19/26 | 24/29 |
| 11369093 | - | WE6/23/7 | 7.5 | 0.10 | 19/26 | 24/29 |
| 11359093 | - | WE6/23/25 | 25.0 | 0.20 | 19/26 | 24/29 |

Receiver for light entrainer, with PTFE GP Rotaflo stopcock

| Catalogue No | Legacy No | Alt. No | Capacity, mL | Graduations, mL | Socket size | Cone size |
|--------------|-----------|------------|--------------|-----------------|-------------|-----------|
| 3 11349093 | - | WE6/23/10T | 10.0 | 0.10 | 19/26 | 24/29 |

Receivers for heavy entrainers

| Catalogue No | Legacy No | Alt. No | Capacity, mL | Graduations, mL | Socket size | Cone size |
|--------------|-----------|-----------|--------------|-----------------|-------------|-----------|
| 12475436 | - | WE7/23/3 | 3.0 | 0.05 | 19/26 | 24/29 |
| 4 12465436 | - | WE7/23/12 | 12.5 | 0.10 | 19/26 | 24/29 |

ACROS
ORGANICS

As a supplier of fine chemicals and organic reagents, Acros Organics continue to improve their product range and services to meet today's requirements of the organic, medicinal, research, analytical or biochemist.

You can find our Structure Search feature on our home page, allowing easy searching of all 34,000+ Acros Organics products.

