

1



Freezers, general purpose, Thermo Scientific Revco® ES and GPS series



For routine sample protection where laboratory space is limited.

ES series

- Digital controller display
- High and low audio and visual temperature alarms
- 16mm access port
- Door locks
- Manual defrost

GPS series

- Digital controller and display
- High and low audio and visual temperature alarms
- 25mm access port
- Casters
- Door locks
- Auto defrost

Technical Specification - Specific

	12938784	12948784	12978784	12988784	12998784	12908794
Type	ES, under bench	ES, upright	ES, freezer/ refrigerator combination	GPS, spark proof, solid door	GPS, large capacity, solid door	GPS, large capacity, solid door
Capacity, L	151	232	109 (freezer) 159 (refrigerator)	400	700	1,400
Temperature range, °C	-10 to -30	-10 to -30	-10 to -30 (freezer) 1 to 11 (refrigerator)	-12 to -25	-12 to -25	-12 to -25
Dimensions, external (w x d x h), mm	595 x 610 x 845	595 x 640 x 1,565	545 x 610 x 1,675	600 x 600 x 1,900	790 x 860 x 1,980	1,400 x 800 x 1,980
Shipping weight, kg	65	90	95	95	130	190

2



Technical Specification - General


Electrical supply 230V, 50/60Hz

	Catalogue No	Legacy No	Alt. No	Description
2	12938784	-	151F-AEW-TS	ES under bench freezer, with 3 baskets
	12948784	-	232F-AEW-TS	ES upright freezer, with 2 shelves and 3 baskets
	12978784	-	263C-AEW-TS	ES freezer/refrigerator combination, with 3 shelves, baskets (refrigerator), 1 shelf, 2 baskets (freezer)
1	12988784	-	F400-SAEW-TS	GPS freezer, spark proof, solid door, with 3 shelves
	12998784	-	F700-SAEW-TS	GPS freezer, spark proof, solid door, with 3 shelves
	12908794	-	F14X-SAEW-TS	GPS freezer, spark proof, solid door, with 6 shelves

Have you thought about.....

Dispensers, for Parafilm®, Fisherbrand®

Stores, dispenses and cuts Parafilm® tape and labels.



Catalogue No	Legacy No	Description	Width, mm	Depth, mm	Height, mm
11937994	-	Acrylic dispenser	119	168	175
11350040	-	ABS dispenser	171	120	144

Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.