

Freezers, laboratory, chest, low temperature



- Store contents down to -45°C
- Fitted lock as standard
- Roller castors for ease of movement
- Internal and external temperature alarm

Available in three sizes, and ideal for use in all environments requiring low temperature storage.

Technical Specification - Specific

	12958674	12968674	12978674
Capacity, L	140	296	383
Dimensions [w x d x h], mm	720 x 600 x 885	1,260 x 600 x 885	1,560 x 600 x 885
Mass, kg	55	76	88

Technical Specification - General

Temperature range, °C	-10 to 45
Refrigerant	R404
Class	Climate Class N
Electrical supply	230V, 50/60Hz

Catalogue No	Legacy No	Alt. No	Capacity, L
12958674	-	CLT140	140
12968674	-	CLT296	296
12978674	-	CLT383	383



Refrigerator/freezer combinations, laboratory



- Spark free
- Fitted lock as standard
- Anti-microbial handle
- Bench top

Technical Specification - Specific

	12978664	12908664	12918664	12988664
Capacity, L	77 (refrigerator) 6.5 (freezer)	106 (refrigerator) 13 (freezer)	159 (refrigerator) 104 (freezer)	190 (refrigerator) 134 (freezer)
Temperature range, °C	2 to 10 (refrigerator) -6 to -9 (freezer)	2 to 10 (refrigerator) -18 to -25 (freezer)	2 to 10 (refrigerator) -18 to -25 (freezer)	2 to 10 (refrigerator) -18 to -25 (freezer)
Refrigerant	R134a	R600a	R600a	R600a
Dimensions [w x d x h], mm	502 x 540 x 660	540 x 600 x 860	540 x 600 x 1,680	600 x 600 x 1,870
Mass, kg	25	34	58	70

Technical Specification - General

Class	Climate Class N
Electrical supply	230V, 50/60Hz

Catalogue No	Legacy No	Alt. No	Description	Capacity, L
12978664	-	ISR27C	Bench top	77 (refrigerator) 6.5 (freezer)
12908664	-	LSC119UK	Under bench	106 (refrigerator) 13 (freezer)
12918664	-	LSC263UK	Freestanding	159 (refrigerator) 104 (freezer)
12988664	-	LSC324UK	Freestanding	190 (refrigerator) 134 (freezer)

