

## Base, slotted



To support flat objects such as boards, manometers, etc., up to 34mm thick in the vertical position and for general laboratory use. Cast metal, finished in acid-resisting paint, with jaws 40mm high and fitted thumb screw. Dimensions of base 95mm x 75mm, mass 1.4kg.




Catalogue No	Legacy No	Description
12328139	-	Slotted base support

## Bases, retort stand, rectangular, Fisherbrand®



Pressed steel finished in acid resistant textured acrylic. 10mm threaded holes. The single hole stand has a hole centrally located on the width. The double hole stand has one hole centrally located on the width and one centrally located along the length.



Catalogue No	Legacy No	Description	Dimensions, mm
11548092	-	Single hole	160 x 100
11558092	-	Single hole	200 x 125
11568092	-	Single hole	250 x 160
11578092	-	Double hole	315 x 200

## Bases, retort stand, rectangular, Fisherbrand®



Bases of blue acrylic finished steel with hole (10mm), threaded at the centre of one short side. Fitted with rubber feet. Without rod.



Catalogue No	Legacy No	Dimensions, mm
11518102	-	160 x 100
11528102	-	200 x 125
11538102	-	250 x 160
11548102	-	315 x 200

For Spoons, please see Scoops pages 1366 to 1368



## Bases, retort stand, rectangular



Bases of white acrylic finished steel with hole threaded at the centre of one short side. Fitted with rubber feet. Without rod.



Catalogue No	Legacy No	Length, mm	Width, mm
11779593	-	160	100
11795376	-	200	125
11705386	-	250	160
11715396	-	315	200

## Stands, burette and funnel

Hardwood, with interchangeable components. Bases are 205mm x 110mm (single) and 290mm x 110mm (double). Rods are 460mm x 16mm (dia). Burette clamps have lined jaws and funnel holders are slotted; all are fitted with a plastic screw and rubber feet on base.



Catalogue No	Legacy No	Description
12906796	-	Burette stand, single
12916796	-	Burette stand, double
12926796	-	Funnel stand, single
12936796	-	Funnel stand, double
12946796	-	Burette/funnel stand

Catalogue No	Legacy No	Description
12956796	-	Burette clamp, single
12966796	-	Funnel support, single
12976796	-	Funnel support, double

## Bases, retort stand, tripod shape



Cast iron tripod spread base painted blue, with rubber feet. The central rod hole is tapped. Without rod.



Catalogue No	Legacy No	Alt. No	Radius, mm
12996796	-	3025/2	100
12338139	-	3025/3	135
12916806	-	3025/4	175