



Magnetic follower sink strainer



A simple solution to prevent magnetic followers being lost down sinks.

- Attracts and holds magnetic followers whilst slits in the side and skirt allows fluids to drain easily
- Flexible, chemical resistant sink strainer with four strong magnets in each corner
- Wide skirt to fit a variety of sink drains
- Four legs on underside anchor strainer into the sink drain

Catalogue No	Legacy No	Alt. No	Colour
11823210	-	F37787-0000	Blue
11833210	-	F37787-0001	Green



Shaft stirrers, PTFE, Fisherbrand®



- Chemically inert
- Easy to clean
- No cross contamination
- Shatterproof
- Non-scratch
- Use in temperatures up to approximately 280°C
- Plus exposed stainless steel core version (EX) for extra rigid clamping in the drive chuck
- Anchor stirrers are for mixing high viscosity substances at slower speeds
- Centrifugal stirrers are for round vessels with narrow necks, the blades open with increasing speed, medium to high speeds are required



PTFE shaft stirrers have a unique, fully moulded construction featuring a rotor moulded to a fully PTFE encapsulated steel reinforced drive shaft.



	Catalogue No	Legacy No	Description	Shaft diameter, mm	Length, mm	Rotor diameter, mm
1	10779372	-	Screw propellor	8	300	40
	10689152	-	Screw propellor	8	400	40
2	10351542	-	Screw propellor (EX)	9.5	650	60
	10279841	-	Screw propellor	10	550	70
	10701493	-	Screw propellor	10	650	70
	10762603	-	Screw propellor (EX)	12	650	80
	10637501	-	Anchor	6	300	80
3	10248020	-	Anchor	6	400	80
	10115590	-	Anchor	6	500	80
	10629542	-	Anchor	8	300	80
	10496861	-	Anchor	8	400	80
	10516712	-	Anchor	9.5	650	100
	10418791	-	Anchor (EX)	9.5	750	140
	10446101	-	Anchor (EX)	10	650	140
	10456101	-	Anchor (EX)	12	650	140
	10084650	-	Centrifugal	6	400	40
	10626351	-	Centrifugal	6	500	50
4	10467441	-	Centrifugal	8	400	40
	10576902	-	Centrifugal	9.5	650	70
	10586902	-	Centrifugal (EX)	9.5	650	70
	10268000	-	Centrifugal (EX)	9.5	750	70
	10457251	-	Centrifugal	10	550	70
	10127801	-	Centrifugal (EX)	12	650	80
	10187701	-	Centrifugal (EX)	12	1,000	80

For Fisherbrand® stirrer guides and shaft guides, see page 1265

