

Thermometers

liquid filled



Thermometers, mercury free, blue spirit filled, Lo-Tox™



Catalogue No	Legacy No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
12389349	-	44/812/8	-10 to 210	1	405
12399349	-	44/809/8	-35 to 50	1	305
12309359	-	44/808/8	-10 to 50	0.5	305
13078583	-	44/810/8	-20 to 110	1	305
13088583	-	44/811/11	-20 to 150	1	305
13038593	-	44/830/7	-10 to 110	1	155
13048593	-	44/850/5	-10 to 50	1	155

Catalogue No	Legacy No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
11907341	-	44/804/8	-10 to 260	1	405
13098583	-	44/801/8	-35 to 50	1	305
13008593	-	44/800/8	-10 to 50	0.5	305
13018593	-	44/802/8	-20 to 110	1	305
13028593	-	44/803/8	-20 to 150	1	305



Thermometers, mercury free, blue spirit filled, Lo-Tox™



General purpose/stirring thermometer with a low toxic, environmentally friendly blue filling, which is easy to read against the white background. Each thermometer is manufactured to traceable standards, and comes with a blue non-roll top and safety suspension ring.

Fractional degree
Total immersion

Catalogue No	Legacy No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
13008603	-	44/821/0	-1 to 51	0.1	405
13018603	-	44/823/0	-1 to 51	0.2	405
12339359	-	44/824/0	-1 to 101	0.2	405

76mm immersion

Catalogue No	Legacy No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
13038603	-	44/820/0	-1 to 51	0.1	405
13048603	-	44/822/0	-1 to 51	0.2	405
13058603	-	44/824/0	-1 to 101	0.2	405



Thermometers, mercury free, brown spirit filling, partial immersion, Lo-Tox™



Mercury free general purpose/stirring thermometers with a low toxic brown spirit filling. This filling enables accurate temperature readings to be taken above 300°C. Each thermometer is manufactured to traceable standards and comes with a blue non-roll top and safety suspension ring.

Catalogue No	Legacy No	Alt. No	Range, °C	Subdivisions, °C	Length, mm
12319359	-	44/815/0	-10 to 300	1	305
13068593	-	44/814/0	-10 to 330	2	305