

### Timer, Clip-It™, Fisherbrand®, Traceable®



- Counts down, alarms, and has stopwatch mode
- Clips to your lab coat, belt, notebook, clipboard, or anywhere
- Six digit LCD display shows the time remaining in hours, minutes, and seconds
- Fast and intuitive five button operation
- Stopwatch mode: 1s to 99h and then repeats, single button operation
- Countdown mode: 1s increments from 100h to 1s
- Alarm: sounds at zero for 60s
- Time-out feature: both up and down counting
- Big digit, 8.5mm high display
- Rugged, chemical resistant ABS plastic case, magnet for placing on metal, and a 1½ year, continuous use, replaceable silver oxide battery
- Size: 51mm x 38mm x 19mm
- Mass: 14g



Catalogue No	Legacy No	Description
11765873	FB70347	Clip-It™ timer

### Timer, memory, Fisherbrand®, Traceable®



- Perfect to dedicate to repetitive tests
- 'Forever memory' returns to last time set with one touch
- Set count down with alarm from 24hr to 1min
- Easy three button operation
- With clip, magnet and stand
- Accuracy: 0.01%
- Dimensions: 60mm (diameter) x 13mm
- Supplied with batteries



Catalogue No	Legacy No	Description
11735873	FB70234	Timer, with memory

#### Spare part

Catalogue No	Legacy No	Description	Pack qty
12326968	BMD-140-011E	Batteries, alkaline, 130mAh, 1.5V (LR44)	4

### Timer, high volume, Fisherbrand®, Traceable®



- Counts down and alarms
- 4 digit LCD 19mm high display shows the time remaining in minutes and seconds
- Countdown: 1s increments from 99m, 59s
- Automatic memory resets display to last programmed time, ideal for repetitive tests
- High volume alarm sounds at zero for 30s
- Magnetic back, flip open easel for standing and a spring fastener for clipping to a lab coat
- Accuracy: 0.01%
- Size: 63mm x 63mm x 13mm
- Mass: 57g
- Supplied with a 2 year, continuous use, AAA battery



Catalogue No	Legacy No	Description
11799795	FB70325	Timer, high volume alarm

### Timer, extra loud alarm, Fisherbrand®, Traceable®



- 100dB alarm may be heard in the noisiest laboratory or plant environments
- Two timing channels may be used simultaneously, countdown/up, timeout function
- Features stopwatch and time of day clock
- Attach with magnetic back or use the flip open stand
- Timing capacity: 100h to 1s
- Accuracy: 0.01%
- Dimensions: 100mm x 63mm x 12.5mm
- Supplied with batteries



Catalogue No	Legacy No	Description
11542793	FB70232	Timer Fisherbrand® with an ISO 17025 A2LA Traceable® NIST certificate

#### Spare part

Catalogue No	Legacy No	Description	Pack qty
12737686	BMT-680-040W	Batteries alkaline, AAA	4

Traceable® products are supplied with a serial numbered certificate, indicating traceability to NIST standards (National Institute of Standards and Technology).